



InterGuide tracheal tube introducers and InterForm intubating stylets

For use in the difficult airway




InterGuide tracheal tube introducer

InterGuide is a flexible tracheal tube introducer commonly known as a bougie. It allows positive location of the trachea and subsequent placement of the endotracheal tube in difficult airway situations. InterGuide is single use, latex free and supplied sterile in easy-to-open individual packs.




▶ Video available at www.intersurgical.com

Atraumatic coude tip

To reduce the potential for patient trauma



Low-friction surface
For easy insertion and operation

Code	Description	Size	Outer Diameter (OD)	Length	Box Qty.
8070015	InterGuide tracheal tube introducer bougie	15FR	5.0mm	700mm	10 
8070010	InterGuide tracheal tube introducer bougie	10FR	3.3mm	700mm	10 
8070006	InterGuide tracheal tube introducer bougie	6FR	2.0mm	530mm	10 

InterForm intubation stylet




InterForm is a malleable stylet, allowing the user to form the endotracheal tube into a suitable shape to ease insertion and give more control of the airway device. InterForm is single use, latex free and supplied sterile in easy-to-open individual packs.

Aluminium core

Allows user to easily form suitable shapes



Low-friction surface
For easy insertion and removal

Code	Description	Size	Outer Diameter (OD)	Length	Box Qty.
8080014	InterForm intubation stylet	14FR	4.7mm	340mm	10 
8080010	InterForm intubation stylet	10FR	3.3mm	340mm	10 
8080006	InterForm intubation stylet	6FR	2.0mm	275mm	10 

 Sterile

IS6.5 US • Issue 3 05.16