

# Trauma Air Pants

## (Anti-shock Trousers)

LSP Trauma Air Pants anti-shock trousers are designed for emergency care for the victim experiencing hypovolemic shock. They provide complete counter-pressure around the legs and abdomen, combining safety and comfort for the patient.

Both the adult and pediatric anti-shock trousers provide the following specialized features:

- Handles are strategically located in the torso area and both legs for easy positioning and patient movement
- All bladders have a 5-year warranty against breakage
- The pant material has a 2-year warranty
- Unique design allows the pants to open outward and lie flat for easy application and patient accessibility
- Color-coded Velcro straps and luminescent pressure gauges for easy nighttime application
- Inflation pump can be used with foot, hand or knee pressure
- Adult pants are easily downsized to fit either a small adult or a pediatric patient
- Soft-pack carry bag is padded for protection, and the strap system converts to a backpack for carrying ease



# Trauma Air Pants (Anti-shock Trousers)



L600

## SPECIFICATIONS

**Outside Fabric:** 420 denier nylon

**Inner Bladder:** Clear 12 mil polyurethane; ultra-violet, bacteria, and fungus stabilized; suitable for extensive water immersion and air retention; 5-year warranty.

**Pump:** High-grade neoprene

**Hoses:** Surgical-quality Latex

**Hose Connectors, Fittings, Adapters:** Injection-molded high-impact moisture-resistant nylon

**Gauges:** Phosphorescent gauges, 0-200 mm Hg

**Pressure Relief Valves:** Pop-off pressure 110-135 mm Hg

**Velcro Closures:** Color-coded white, black

**Bag:** 200 denier nylon soft side

## ORDERING INFORMATION

- |               |  |
|---------------|--|
| <b>L600</b>   | Trauma Air Pants, adult, with 3 gauges, pump assembly, bag                   |
| <b>L600-G</b> | Trauma Air pants, adult, with 3 gauges, pump assembly, bag (Germany)         |
| <b>L601</b>   | Trauma Air Pants, adult, with 1 gauge, pump assembly, bag                    |
| <b>L602</b>   | Trauma Air Pants, adult, with pressure relief valves, pump assembly, bag     |
| <b>L610</b>   | Trauma Air Pants, pediatric, with 3 gauges, pump assembly, bag               |
| <b>L611</b>   | Trauma air Pants, pediatric, with 1 gauge, pump assembly, bag                |
| <b>L612</b>   | Trauma Air Pants, pediatric, with pressure relief valves, pump assembly, bag |



Authorized EU Representative  
Medizintechnik & Consulting  
Effengrube 14  
23552 Lubeck, Germany  
TEL: +49-451-8090-4200  
FAX: +49-451-8090-4231



Allied Healthcare Products, Inc. ▪ 1720 Sublette Avenue ▪ St. Louis, MO 63110 USA  
314/771-2400 ▪ 800/444-3954 ▪ FAX: 314/771-0650 ▪ www.alliedhpi.com

U.S.  
Toll Free: 800/444-3940  
FAX: 800/477-7701

Canada & Mexico  
Toll Free: 800/446-0552  
FAX: 800/246-6201

International  
Telephone: 314/268-1683  
FAX: 314/771-5183

ISO 13485: 2003

©Allied Healthcare Products, Inc.®

All specifications are nominal and  
subject to change without notice.

Form 90-00-0020  
December 2011