



Your Need ... Our Innovation[®]

ONLY ONE

naso-flo[®]

Nasopharyngeal Airway Device

- Oxygen port and elbow allow for direct connection to O₂ source (higher FiO₂ % delivery)
- Soft silicone material designed for less traumatic insertions
- Humidification vents allow for heat/moisture transfer on inhalation
- Color-coded elbows aid in correct size selection
- Optional filtered respiratory indicator (CO₂ monitor connection)

naso-flo w/Oxygen Port 10/Box

Item #	Inside Diameter	French Size	Outside Diameter	Length
JQ-NF104	4.0mm	18FR	6.0mm	100mm
JQ-NF105	5.0mm	22FR	7.3mm	115mm
JQ-NF106	6.0mm	25FR	8.3mm	130mm
JQ-NF107	7.0mm	29FR	9.6mm	140mm
JQ-NF108	8.0mm	32FR	10.7mm	170mm
JQ-NF109	9.0mm	36FR	12.0mm	190mm

naso-flo w/Oxygen Port/EtCO₂ Connector w/Filter 10/Box

Item #	Inside Diameter	French Size	Outside Diameter	Length
JQ-NF204	4.0mm	18FR	6.0mm	100mm
JQ-NF205	5.0mm	22FR	7.3mm	115mm
JQ-NF206	6.0mm	25FR	8.3mm	130mm
JQ-NF207	7.0mm	29FR	9.6mm	140mm
JQ-NF208	8.0mm	32FR	10.7mm	170mm
JQ-NF209	9.0mm	36FR	12.0mm	190mm

Patent Pending

6 Sizes Available



naso-flo with Oxygen Port

naso-flo with Oxygen Port/EtCO₂ Connector with Filter

11300 - 49th Street North • Clearwater, Florida 33762-4807 USA

Tel. (800) 835.6633 • (727) 573.0088 • Fax (800) 990.6375

www.mercurymed.com

