



THE EMERGENCY BANDAGE®

- Stop n' Go Release
- FDA Approved
- Internationally Patented
- Specified for Use by the US Army
- Enhanced Vacuum-Sealed Packaging

MILITARY COMMERCIAL & INDUSTRIAL SITES FIRST AID KITS FIRST RESPONDERS HOMELAND SECURITY LAW ENFORCEMENT DISASTER PREPAREDNESS EMS



www.ps-med.com
info@ps-med.com

5310 Elm Street
Houston, Texas 77081

Tel 713-723-6000
Toll Free 1-888-737-7978



MOBILE PADS **STOP N' GO** PRESSURE APPLICATOR **CLOSURE BAR**



The Emergency Bandage®

The combat proven bandage of choice for U.S. and international military units.

This life saving bandage is the ideal solution for emergency treatment when every second counts.

Its innovative design staunches blood flow for hemorrhagic wounds in emergency situations.

Features and Advantages


Stop n' Go Release



Vacuum Sealed



Pressure Applicator



Mobile Pad



All-In-One Device consolidating a primary dressing, pressure applicator, secondary dressing, and a foolproof closure bar into one compact unit

New Stop n' Go Release prevents unrolling and ensures wrapping control

Sterile, non-adherent pad and pressure applicator applies pressure to any site

Allows an injured person to self-apply the bandage with one hand

Flexible and versatile for nearly all traumatic extremity injuries

Product Numbers and Sizes:

NSN# 6510-01-460-0849	part# FCP-01	4" width
NSN# 6510-01-492-2275	part# FCP-02	6" width
NSN# 6510-01-515-7528	part# FCP-05	6" width + mobile pad
NSN# 6510-01-532-6656	part# FCP-09	12" x 12" abdominal pad
	part# FCP-07	4" width + mobile pad



www.ps-med.com
info@ps-med.com

