## Zenith ALS/BLS or Airway Bag

The Zenith ALS/BLS or Airway bag is the perfect blend of size and space. It is ideal for the first in bag, featuring both backpack straps and over the shoulder option.

Made to accommodate a "D" sized tank and multiple sizes of BVM's, Intubation, IO, Pulse Oximetry, and other initial assessment and wound treatment supplies.



## Features:

- 3 exterior pockets, and 2 glove stuff pockets
- Large interior compartment
- Removable sectional dividers for keeping equipment organized
- 2 large removable clear vinyl fronted pockets for storage and organization
- Padded base for protected storage of 02 tank
- Nylon straps and side release buckles to ensure secure storage

## **Superior construction:**

- Made with reversed zippers to avoid bending zipper teeth
- Zippers are located not on the edges of the bag, so that they don't rub when corners and edges are beat up from storage, and dragging on the floors of ambulances and Fire trucks
- Only #10 YKK zippers the best in the industry
- Made with 2 handle top carry system to avoid zipper stress when hefting or carrying the bag with the hand handles.
- All 1000 D Cordura
- High Visibility Zipper pulls for all exterior compartments and pockets
- Tuck-Away back pack straps for hands free carrying

New Thomas EMS Fluid free material wears like Cordura, keeps fluids from penetrating material.

ltem	Description	Price
TT1280	Zenith ALS/BLS or Airway Bag, Red	\$350.00
TTFF1280	Zenith ALS/BLS or Airway Bag, Black, Fluid Free Material	\$450.00





Thomas EMS Toll Free (800) 445-3640 (801) 262-6503 • Fax (801) 268-9272 www.thomasems.com



