

# Ambu<sup>®</sup> Man School



## Key features

- Ambu hygienic system to prevent cross contamination
- Compact and portable training manikin
- Easy training of chest compressions
- Closed torso
- Easy placement of AED pads
- 10 manikins for optimal class training

## AmbuMan School

AmbuMan School is a training manikin which offers human appearance. It is dedicated to the training of all procedures for basic life support resuscitation.

The manikin can be used for realistic training of external chest compressions and the placement of AED pads in order to get comfortable with AEDs.

The unique hygienic system prevents cross contamination when demonstrating and training the principals of ventilation through mouth and nose.

The AmbuMan School is intended to be used for education of pupils in schools under the supervision of an instructor/teacher.

### Conforms to Guidelines

The AmbuMan School allows instructors to perform CPR training following the latest CPR guidelines.

### Compact and portable

Lightweight and compact, the AmbuMan School is easy to transport and set up.

### Easy AED training

Due to the closed shaped torso, defibrillation pads can be placed and simulated AED training can be performed.

# Specifications

Weight	0.8 kg
Length	45 cm
Max. compression depth	6 cm

# Material

Product	Material
Straps	Nylon
Base plate	Polycarbonate
Bellow	Thermoplastic Elastomer
Chest	Polypropylene
Head	Silicone
Face piece	Polyvinylchloride
Lung bag	Polyvinylchloride
head bag	Polyethylene

# Order information

Item no.	Description
A268 007 000	AmbuMan School (Set of 10)
268 000 571	Face piece AmbuMan School (Set of 10)
268 000 572	Head bag AmbuMan School (Set of 100)
268 000 573	Accessories AmbuMan School (Set of 10)



Closed torso



Placement of AED pads



Ambu hygienic system



Accessories

**Ambu Inc.**  
6230 Old Dobbin Lane  
Columbia, MD 21045  
Tel. 800 262 8462  
Fax 800 262 8673  
ambuUSA.com

**CE** US: Rx only