

# O2-RESQ™


## 3-SET™ O2-CPAP™ Valve


The Pulmodyne® 3-SET O2-CPAP Valve is a new cutting-edge device that allows the clinician to titrate the level of CPAP with the simple turn of a dial. Eliminating the need to carry extra equipment or supplies.

**Simple**

**Effective**

**Titration**

- *Pre-assembled O2-RESQ System includes:*
  - Flow Generator with Filter
  - Expandable Circuit
  - 3-SET O2-CPAP Valve
  - Mask and Head Strap
- *BiTrac ED Full Face Mask features an OmniClip™ with multi-positioning Silicone Forehead Pad*
- *Flow Generator provides 30% O<sub>2</sub> concentration at a Continuous Flow up to 140 l/m*
- *3-SET O2-CPAP Valve with adjustable pressure settings from 5.0, 7.5, 10.0cm H<sub>2</sub>O*
- *3-SET O2-CPAP Valve provides fixed pressure at flow rates from 60 l/m*
- *Easy to use adjustable valve*
- *Eliminates need for a manometer*
- *Control over protocols*
- *Cost Effective*
- **DISPOSABLE** - Single Patient Use
-  **Latex Free**



*Complete line of pre-assembled disposable CPAP Systems for Immediate/Emergency Ventilation*

# O2-RESQ™

## 3-SET™ O2-CPAP™ Valve

Part No.	Product Description	Size	Qty/Case
<b>Systems</b>			
<b>O2-RESQ System: BiTrac ED Full Face Mask with 30/22mm Elbow, OmniClip and Head Strap, Flow Generator and Expandable* Anti-Asphyxia Circuit with CPAP Valve and Adapter (22mmM x 30mmM)</b>			
313-7055X	BiTrac ED Mask with 3-SET CPAP Valve (5.0, 7.5, 10.0cm)	Adult Large	10
313-7056X	BiTrac ED Mask with 3-SET CPAP Valve (5.0, 7.5, 10.0cm)	Adult Medium	10
<b>*Expandable Anti-Asphyxia Circuit expands to 72"</b>			
<b>**Contact Pulmodyne for Additional O2-RESQ System Configurations available.</b>			
<b>Accessories</b>			
<b>3-SET O2-CPAP Valves with Connectors (22mm and 30mm)</b>			
313-7725	3-SET CPAP Valve (5.0, 7.5, 10.0cm)		10
<b>Valved "T" Adapter with one arm 22mmM and one arm 30mmM</b>			
313-2230			30
<b>O2-RESQ Filter (30mmM x 30mmF; BFE and VFE 99.9999%)</b>			
313-7043			10



BiTrac ED Mask



3-SET O2-CPAP Valve



O2-RESQ Filter



Valved "T" Adapter

### Masks

#### BiTrac ED™ Full Face Mask and Head Strap\*

313-7031	Adult Extra Large
313-7030	Adult Large
313-7029	Adult Medium
313-7028	Adult Small

\*All Mask are 10 per case



... bringing change to life®

Pulmodyne  
2055 Executive Drive  
Indianapolis, IN 46241 USA  
T: 317.246.5505  
F: 317.246.5501  
www.Pulmodyne.com

Pulmodyne and *bringing change to life* are registered trademarks of Pulmodyne, Inc. O2-RESQ, BiTrac ED, BiTrac Shield ED, BiTrac MaxShield ED, OmniClip, O2-CPAP and 3-SET are trademarks of Pulmodyne, Inc.