

▲ COBALT®

POWDER-FREE EXAMINATION GLOVES

BLUE **NITRILE NON-STERILE AMBIDEXTROUS FULLY TEXTURED BEADED CUFF** FOR SINGLE USE ONLY **NON-LATEX**

This soft and stretchy nitrile glove is the best-selling glove from Microflex® brand of products. Durable and puncture resistant witha fully textured surface for enhanced grip. Ideal for jobs requiring exceptional tactile sensitivity.

- Unique, cobalt blue coloring
- Consistent product with exceptional fit and feel





PRODUCT CODE				
XS	N190			
S	N191			
М	N192			
L	N193			
XL	N194			

PACKAGING			
100	Gloves/Box		
10	Boxes/Case		

GLOVE PROFILE		ASTM D 6319	COBALT®
Length	mm (inches)	min. 220 (8.7)	240 (9.5)
Cuff Thickness	mm (mils)	n/a	0.09 (3.5)
Palm Thickness	mm (mils)	min. 0.05 (2.0)	0.10 (3.9)
Finger Thickness	mm (mils)	min. 0.05 (2.0)	0.12 (4.7)
Tensile Strength	Before Aging After Aging	min. 14 min. 14	18 14
Elasticity/Elongation	Before Aging After Aging	min. 500 min. 400	500 400
Pinhole Defect Rate*	(AQL)	max. 2.5	1.5

^{*} Pinhole Defect Rate is FDA minimum requirement.

















