## V.E.R.S.A. PROINFECTION CONTROL COMPLETE

Versatile Emergency Response System Assist



M4L Ballistic Armored

The V.E.R.S.A.™ PRO X is the most versatile emergency response pack ever created with COMPLETE INFECTION **CONTROL.** Constructed internally and externally from non-absorbing wipe-clean materials, the V.E.R.S.A.™ PRO X is specifically designed to easily clean and remove blood borne pathogens while offering a completely customizable backpack capable of securing internal or external full size or mini modules, is MOLLE compatible, has adjustable internal dividers, can carry an M9 oxygen cylinder, and can include M4L Ballistic Armored Protection. The V.E.R.S.A.™ PRO X delivers a highly versatile, organized solution for your equipment so you can respond faster and safer.

Low profile COMFORT-GRIP top carry handle with black logo

Oversized YKK zippers

TS2 Ready<sup>™</sup> hideaway zippers for 2 optional side TS2 Ready™ modules



M4L Ballistic Armor provides optional 12" (tall) x 10" (wide) NIJ Level IIIA armor protection to vital chest region





Fully concealable waist and backpack straps when not in use

Durable bottom skid runners



Drain hole bottom

Loop tape lined panel attaches a variety of internal modules



Supports up to M9 Oxygen Cylinder

Main pocket opens complete



Zippered clear document pocket

TS2 Ready™

panels

hideaway zippers



Front compartment

supplies and small gear

organizes office

Large internal compartment can be customized with dividers or internal modules




Infection control materials designed to easily clean and remove blood borne pathogens

Black matte finish metal hardware

Padded, vented backpack system with contoured yoke harness system easily distributes weight for maximum comfort

3" x 11" front name tape area, MERET tape included



Front MOLLE panel for accessory customization

Zippered internal mesh compartments



Overlapping internal panel secures contents



Adjustable sternum and waist belt



Paneled with laser cut MOLLE to support a TS2 Ready™ module or any MOLLE compatible gear

Product Dimensions: 18.5" x 13" x 10.5" Product Weight: 6lbs, 15oz standard M4L Lifetime Warranty

## \*M4L BALLISTIC CAPABILITIES

NIJ Standard 0108.01 Threat Protection (Level IIIA)

- 9mm Luger 124gr FMJ
- .44 MAG 240gr Lead SWC GC

Special Threat Protection

- Federal 9mm Luger 124gr Hydra-Shok JHP
- Federal .40 S&W 165gr Hydra-Shok JHP
- Federal .45 AUTO +P 200gr GDHP

Ballistic Technology Designed for MERET by Angel Armor



Modules and supplies shown in the images above are not included.