

PRECISION MEDICAL®

Oxygen Regulators



EASY DIAL REG

OXYGEN REGULATOR

The world's most compact and convenient regulators!

They're Compact and Lightweight...

The Easy Dial Reg is only 4-1/4" long and weighs a mere 7.8 oz. as a result of the offset gauge and new multi-spring technology. The Easy Gauge Reg is only 4.63" long and weighs 12 oz. with the contents gauge on the front and the flow gauge on the top.

They're User-Friendly...

Both types of regulators have easily accessible flow adjustment knobs making flow adjustment a cinch.

They're Safety is Built in...

The Easy Reg line of oxygen regulators feature a built-in pressure relief valve as well as a protected contents gauge. The low force Viton yoke seal provides a better tank-to-reg seal.

They're Safe...

The Easy Dial Reg has passed the **ASTM G175 Ignition Sensitivity and Faulty Tolerance Test.**

They're Accurate...

Unique design minimizes normal flow variation that occurs as tank pressure drops.


They're Flexible & Versatile...

The Easy Dial Reg comes in two (2) body designs, yoke or nut style, and is available in five (5) flow ranges: 31-4000 cc/min, 1/32-4 lpm, .25-8 lpm, .50-15 lpm and .25-25 lpm. The Easy Gauge Reg also comes in the two (2) body designs and two (2) flow ranges: 0-8 lpm and 0-15 lpm.



EASY DIAL REG


OXYGEN REGULATOR



168715DSDS shown

EASY GAUGE REG

OXYGEN REGULATOR



168715GE shown

SPECIFICATIONS (EASY DIAL REG)

Flow Range	Position Increments	Yoke Style	Nut Style
		(870)	(540)
0-4000 cc/min	OFF, 31, 62, 125, 250, 500, 750, 1000, 1500, 2000, 3000 and 4000	168704DCL	165404DCL
0-4 lpm	OFF, 1/32, 1/16, 1/8, 1/4, 1/2, 3/4, 1, 1-1/2, 2, 3 and 4	168704DL	165404DL
0-8 lpm	OFF, .25, .50, 1, 1.5, 2, 2.5, 3, 4, 5, 6 and 8	168708D	165408D
0-15 lpm	OFF, .50, 1, 1.5, 2, 3, 4, 5, 6, 8, 10 and 15	168715D	165415D
0-25 lpm	OFF, .25, .50, 1, 2, 3, 4, 6, 8, 10, 15 and 25	168725D	165425D

ORDERING INFORMATION

Easy Dial Reg

DSDS	DSDR	DSR	NG	0706
168715DSDS shown Outlet DISS O2	168715DSDR shown Outlet DISS O2 90° pointing down	168715DSR shown Outlet Hose Barb 90° pointing down	168708DNG shown No Gauge	168715DSDS/0706 shown Extension used with DSDS for humidifier bottles

Easy Gauge Reg

Flow Rate	Yoke Style (870)	Nut Style (540)
0-15 lpm	168715G	165415G
	168715GE	165415GE



ISO 13485
Certified

Precision Medical, Inc.
300 Held Drive • Northampton, PA 18067
Phone: 610-262-6090 Toll-Free Phone: 1-800-272-7285
Fax: 610-262-6080 Toll-Free Fax: 1-800-352-1240
Visit us on the Web: www.precisionmedical.com