

Thomas™ Tube Holder


Tho



All Laerdal Thomas Tube Holders share features that set them apart from other airway tube holders. For instance, the Laerdal Thomas Tube Holder is the only tube holder on the market that has a quick-set screw for a secure hold after tube placement. Tape is not necessary.

Laerdal Thomas Tube Holder is a unique (and better) tube holder over others on the market, no matter what size or color you choose.

- Eliminates the need for cumbersome taping or other fixatives
- Slide stick facilitates fast, efficient application of the head strap around the patient's neck without movement of the head
- Large opening allows for application of a wide range of single and double-lumen airway tubes
- Ventilations may continue during application, adjustment and removal of the Thomas Tube Holder



HEALTHCARE



EMS



MILITARY

mas Tube Holder



Thomas Tube Holder Features

Wide Padded Velcro® Head Strap Features:

- Provides secure and comfortable fit, even when wet
- Color-coded for easy identification between adult and pediatric
- Will not stretch

Military

Adult

Pediatric

Suctioning Features:

- Larger access opening to allow for ease of suctioning and visual inspection of mouth

Slide Stick Features:

- Facilitates fast, efficient application of the head strap around the patient's neck without movement of the head

Security Features:

- Set Screw mechanism quickly secures ET tubes and allows for fast readjustment

Enlarged Access Port Features:

- Accommodates multiple size ET tubes, as well as Combitubes®, King Tubes, LMAs, and LMA Supreme™
- BVM does not need to be removed to attach the Thomas Tube Holder
- Allows for in-use suctioning of the oropharynx

