

Seeing is Believing




INTRODUCING

CO₂ Easy™

CARBON DIOXIDE DETECTOR

Optimize Visibility, Clarity and Confidence in CO₂ Detection

Featuring the largest visual indicator available, the CO₂ Easy™ is ideally mounted up front where clinicians and first responders can effectively view critical carbon dioxide levels without looking away from the patient. The nontoxic, pH sensitive indicator uses familiar purpleto-yellow color conversion when CO₂ is exhaled.



Part Nur	nber Description	Case Quantity
562134	CO ₂ Easy Carbon Dioxide Dete	ctor10
562135	Pedi CO2 Easy Carbon Dioxide	Detector 10
562136	BagEasy w/ Adult Mask and C	
562233	BagEasy w/ Peep, Adult Masl	c and CO2 Easy6

- CO₂ Easy is a trademark and Bag Easy is a registered trademark of Westmed, Inc.
- © 2006 Westmed, Inc. All rights reserved. PN 74558/Rev. 04

- Built in bacterial filter protects both patient and clinitian.
- C02 Easy for patients weighing 15kg (33lbs.) or above.
- Pedi C02 Easy for patients weighing 1 15kg (2.2 - 33lbs.).
- 3 year shelf life.
- Designed for fast, easy airway assessment.
- Verifies proper endotracheal tube placement.
- Helps prevent unwarranted intubations.
- Confirms cardiac output during CPR.
- Expansive 2" Diameter Visual Indicator makes rapid breath-to-breath response easy to see.
- Anti-fogging agent reduces susceptibility to moisture and humidity, while ensuring continuous lens clarity. (C02 Easy only)
- Complies with Standard of Care Guidelines from the American Heart Association, American College of Emergency Physicians, American Association for Respiratory Care and the American Society of Anesthesiologists.





Westmed, Inc.

5580 S. Nogales Hwy. Tucson, AZ 85706 Phone: 800-975-7987 Fax: 520-294-6061 www.westmedinc.com