



MERET



MERET


877.222.0200
www.MeretUSA.com
info@meretusa.com
Ph: 714.259.4700
Fx: 714.259.4749

MERET OXYGEN REGULATORS

MERET oxygen regulators incorporate the latest in technology to deliver the highest level of safety, reliability, durability, and precision flow. EMS series unibody regulators are constructed from high strength solid brass with brass internal high pressure flow components to ensure added strength and safe operation.

MEDICAL * EMERGENCY * EQUIPMENT & TECHNOLOGY

MAGNUM SERIES ALL-BRASS OXYGEN REGULATOR



Product Dimensions:
6.5" x 3.5" x 3.5"
 Product Weight:
2 lbs 1 oz

MERET CGA REGULATOR OPTIONS

REGULATOR	MODEL	FLOWS (L/MIN)	50 PSI POWER TAKE-OFFS (PTO)	CONSTRUCTION	WARRANTY
MAGNUM	EMSREG8725-xxx	0-25	2 + Barb	Solid 360 Brass, Chrome Plated	LIFETIME
BRUUT	EMSREG8725-xxx	0-25	0-2 + Barb	Solid 360 Brass, Chrome Plated	5 yr
ECONOMY	AREG87xx-xxx	0-4, 0-8, 0-15, 0-25	0, 2 + Barb	Aluminum/ Brass Internal, Anodized	5 yr
MINI	AREG87xxME-BU	0-8, 0-15	Barb Only	Aluminum/ Brass Internal, Anodized	5 yr
FLOWMETER	MFMC-25	0-25	N/A	Aluminum/ Brass Internal, Anodized	5 yr
PRESSURE REDUCER	PRREG54-2D	50 PSI Preset	2	Solid 360 Brass, Chrome Plated	5 yr
OMNIVALVE	AI0-3-xx	0-8, 0-5, 0-25, integral valve	1 + Barb	Solid 360 Brass, Chrome Plated	5 yr

All regulators are individually calibrated and tested 100% for flow and leak. ASTM G175, ISO 13085, ISO 9001:2000

BRUUT

- All brass unibody, metal flow knob
- 0-2 PTO and barb

ECONOMY

- Aluminum outer body
- Multiple flow options + PTO

MINI

- Smallest O₂ reg at only 4.0" in length, 0.7 LB + PTO

FLOWMETER

- Compact, lightweight
- Uses 1/8" Male NPT connection

PRESSURE REDUCER

- All brass, compact
- 2 PTO w/ check valves

OMNIVALVE

- Valve integrated regulator
- All brass, 1 PTO + barb

MEDICAL * EMERGENCY * RESPONSE * EQUIPMENT & TECHNOLOGY