

Everywhere Your First-In Bag Goes...



...The Quickdraw Can Be There

The SSCOR Quickdraw easily fits into an airway or first-in bag so you can have quality medical suction available at the patient's side. The SSCOR Quickdraw provides suction as efficient as a large aspirator but weighs only 2.6 pounds (1.18kg).



Powered by alkaline battery technology!

Alkaline battery technology is ideal for small portable devices carried in kits. Alkaline batteries have a long shelf life, require no charging, pack a lot of power for long run times and are widely available and inexpensive.

SSCOR offers an alkaline battery pack which will power the Quickdraw for 3 hours of continuous use. And now SSCOR also offers the all new 'AAA' battery holder which accepts standard off the shelf 'AAA' alkaline batteries.

SSCOR

The Quickdraw is the ideal aspirator to carry in your 'first in' kit and have on hand every time you are with the patient.

Quickdraw! Have it ready to use at your side on every call!



Items Included



Canister with Hose Barb

The Quickdraw canister with hose barb will accept standard suction tubing. For use with common suction tips including tracheal catheters.

#2488 Each

#2488-20 Case of 20



SSCOR HI-D "Big Stick" Suction Tip

A large bore suction tip (.25" I.D.) with thumb port, allowing the patients airway to be cleared more effectively.

#44241 Each

#44241C Case of 50



Battery Pack

Run time for the SSCOR alkaline battery pack is three hours +/- 10%.

#80613-100 Each



Connecting Tubing 3 Feet

3' suction connecting tubing.

#43203 Each

#120-00001 Case of 50



Battery Holder

Accepts 10 'AAA' batteries. Run times with a SSCOR battery holder will vary from 60 to 100 minutes depending on the brand and freshness of batteries used. (Batteries not included.)

#80615-100 Each

Optional Products



SSCOR DuCanto Suction Catheter

A large bore suction catheter (.26" I.D.) designed for use during routine and emergency airway management.

#200-00002 Each

#200-00002C Case of 40

Specifications

Size:	10.5" L x 4.5" H x 4.35" W (27cm L x 11cm H x 11cm)
Weight:	2.6 lbs or 1.18 kg
Performance:	Maximum negative pressure >500mmHg; Low negative pressure setting approximately 80 - 100mmHg
Collection Canister:	#2488 300 cc/ml capacity; 99.98% bacterial effective filter
Run Time:	#80613-100 battery pack: 3 hours +/- 10% #80615-100 battery holder: 60 - 100 minutes with 10 new 'AAA' batteries (batteries not included)

Visit shop.sscor.com for replacement parts, disposables and accessories.

Complete product information can be found at www.sscor.com.

