

Microflex®
ULTRA ONE®
 Powder-Free Latex
 Extended Cuff Exam Gloves

Item #
 Small UL-315-S
 Medium..... UL-315-M
 Large..... UL-315-L
 Extra Large..... UL-315-XL



Recommended for use in the following industries:
 Automotive / EMS / Industrial Safety / Laboratory

For the highest level of protection and durability in latex examination gloves for high-risk environments, rely on **Ultra One®** gloves. Strong, durable and comfortable, **Ultra One®** gloves offer superior defense against many hazards that you come into contact with daily. With an extended cuff for increased protection, the **Ultra One®** glove is a proven performer.

Ultra One® gloves are the best choice in a latex for heavy-duty/high-risk environments.

- Recommended for use when:**
- Added protection is required for extreme environments
 - Extra coverage over the wrist and forearm is needed for enhanced safety and security
 - Strength is required for your toughest tasks

NFPA certification on file.

Microflex® Ultra One® Profile			For current specification information, please visit www.microflex.com .	
	ASTM D 3578	Ultra One®	Material: Latex	Type: Non-sterile
Length (mm/in):	230 / 9.1	300 / 11.8	Size: S - XL	Color: Natural
Cuff Thickness (mils):	n/a	6.7	Shape: Ambidextrous	Usage: For single use only
Palm Thickness (mils):	min. 3.1	9.8	Cuff: Extended Cuff length (12"), Beaded	Glove Exterior: Textured Fingers
Finger Thickness (mils):	min. 3.1	11.8	Glove Interior: Powder-Free, Chlorinated	Packaging: 500 gloves/case, by weight; 10 boxes/case
Tensile Strength:			50 gloves per dispenser box, by weight	
Before Aging (MPa):	min. 18	33	Caution: This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.	
After Aging (MPa):	min. 14	25		
Elasticity/Elongation:				
Before Aging (%):	min. 650	750		
After Aging (%):	min. 500	600		

The reason is the quality... the difference is Microflex®.
 See for yourself why Microflex® gloves are the first choice of exam glove users nationwide. For **FREE SAMPLES**, please visit www.microflex.com. For the name of a distributor in your area, call: **(800) 876-6866**.

