

Material Safety Data Sheet

PDI® - PROFESSIONAL DISPOSABLES INTERNATIONAL
THE HEALTHCARE DIVISION OF NICE-PAK PRODUCTS, INC
TWO NICE-PAK PARK
ORANGEBURG NY 10962-1376
845-365-1700

REVISION DATE: 3/03

I. PRODUCT IDENTIFICATION

PRODUCT/TRADE NAME: **SANI-DEX® ALC Antimicrobial Alcohol Hand Wipes – Item# Q20372**

HAZARD RATING(NFPA)

HEALTH: 1
FLAMMABILITY: 3
REACTIVITY: 0
SPECIFIC: NONE

EMERGENCY OR INFORMATION TELEPHONE NO: 845-365-1700 (M-F daytime) at other times, contact the local Poison Control Center.

CHEMICAL NAME: Proprietary Mixture

II HAZARDOUS INGREDIENTS PER 29 CFR 1910.1200

Hazardous Ingredients	% w/w	ACGIH TLV	CAS NUMBER
SD ALCOHOL 40 (190 Proof)	65.3	1000 PPM	64-17-5

III-PHYSICAL/CHEMICAL CHARACTERISTICS

COLOR/ODOR/APPEARANCE: Towellette saturated with a clear, colorless liquid with an alcohol odor

BOILING POINT: N/A

FLASH POINT: 74 °F

VAPOR DENSITY: N/A

EVAPORATION RATE: N/A

SOLUBILITY IN WATER: Complete

SPECIFIC GRAVITY (H₂O = 1): 0.890

IV FIRE & EXPLOSION HAZARD DATA

FLASH POINT (Method Used): 74 °F Tag Closed Cup : LEL: N/A UEL: N/A

EXTINGUISHING MEDIA: Carbon dioxide or alcohol type foam

SPECIAL FIRE FIGHTING PROCEDURES: Handle as flammable liquid. Use respiratory protection for fire fighting personnel.

UNUSUAL FIRE AND EXPLOSION HAZARDS: NFPA Class 3 Flammability

V-REACTIVITY DATA

STABILITY: Stable

CONDITIONS TO AVOID: None Known

INCOMPATIBILITY: None

HAZARDOUS DECOMPOSITION OR BYPRODUCT: None Known

POLYMERIZATION: Will not occur

CONDITIONS TO AVOID: None known

SANI-DEX® ALC Item #Q20372

VI-HEALTH HAZARD DATA

EFFECTS OF OVER EXPOSURE:

SKIN	Formulated for topical application
EYES :	Splashes may cause irritation
INHALATION:	None known
INGESTION	Not a normal route of entry

EMERGENCY AND FIRST AID PROCEDURES:

SKIN CONTACT:	Discontinue use if rash
EYE CONTACT:	Flush with cool water for 15 minutes
INHALATION:	None
INGESTION:	Consult Physician
TARGET ORGANS:	None known

VII-SPILL AND DISPOSAL PROCEDURE

SPILL CONTROL:	Pick up with absorbent media
WASTE DISPOSAL METHOD:	Obey local, state and federal regulations
HANDLING AND STORAGE:	Store away from ignition sources such as heat, sparks, open flame, etc...

VIII-CONTROL MEASURES/PROTECTION

RESPIRATION:	None required
VENTILATION:	Recommended
PROTECTIVE GLOVES:	None required
EYE PROTECTION:	Only if splashing is expected
HYGIENIC PRACTICES:	Good housekeeping practices should be followed
OTHER:	None

IX TRANSPORT/SHIPPING

DOT SHIPPING NAME:	Not applicable
TECHNICAL SHIPPING NAME:	Not applicable
DOT SHIPPING CLASSIFICATION:	CONSUMER COMMODITY
DOT ID NO.:	Not applicable
DOT LABEL REQUIREMENTS:	ORM-D
UN/NA NUMBER REGULATIONS:	Not applicable
REPORTABLE QUANTITY:	Not applicable

X DISCLAIMER

THE INFORMATION FURNISHED HEREIN IS BELIEVED TO BE ACCURATE AND REPRESENTS THE BEST DATA CURRENTLY AVAILABLE TO US. NO WARRANTY, EXPRESSED OR IMPLIED IS MADE AND NICE-PAK PRODUCTS, INC. ASSUMES NO LEGAL RESPONSIBILITY OR LIABILITY RESULTING FROM ITS USE.

N/A = NOT AVAILABLE