



Stifneck Select Extrication Collars

Laerdal Spinal Immobilization Solutions

All the adult Stifneck sizes in one easily adjustable, fully assembled collar. Just measure the patient, select the size, adjust, and lock – four simple steps!

The Select collar fits all adult patients and delivers the same superior immobilization you expect from Stifneck collars and, for the smallest patients, the Pedi-Select collar offers an anatomically correct collar for infants and smaller children.



Features

Locking Mechanism: • Four positions: No-Neck, Short, Regular, and Tall • Clearly marked and easy to identify · Ensures stability of collar Audible response for safety • Tactile response for safety

- · Assures symmetrical adjustment



- Molded-in instructions on the back panel
- Unique, large rear panel allows for:
- c-spine palpation
- drainage
- visualization while promoting ventilation for better patient comfort



Tracheal Opening:

- · Oversized tracheal opening allows for easy access for advanced airway procedures
- · Optimal visualization for assessment
- · Easy carotid pulse monitoring

Integrated Cannula Hooks:

- Two positions: left and right
- Properly secures supplemental oxygen devices for fuss-free transportation

Pedi-Select Features: Easily identifiable pink strap

· 3 positions: pedi, midi, mini

980010 Stifneck Select Extrication Collar

980012 Stifneck Select Extrication Collar, Green

980020 Stifneck Pedi-Select Extrication Collar

Benefits

- Four adult collars in one/three pediatric collars in one saves room in bags and ambulance
- Easy-to-use sizing measurements (finger technique)
- Locking mechanism ensures stability of collar and symmetrical adjustment
- Directions on collar make training easy and guarantees instructions always 'on hand' for rescuer
- In-service training video minimizes training time and cost

- Flat storage requires minimum space
- Preassembled assures speedy application
- Oversized tracheal opening allows optimal visualization for assessment and treatment
- Constructed from latex-free quality materials
- Radiolucent and MRI compatible



HOSPITAL



F M S



MILITARY

For additional information or to order, contact your Local Authorized Laerdal Distributor!

