



A BETTER ELECTRODE IS ALL WE DO.



Features and Benefits of the Heart Sync Product Line

Cost Savings

On average we can save an organization up to 15% on their annual defib pad expenditures.

Standardization

We can standardize an organization to one pad from one vendor no matter what make or model of defibrillator you have.

Three Year Shelf Life

Reduce cost and waste. Our Adult and Pediatric pads in low use areas will provide an additional 50% shelf life versus the rest of the market.

Longest Leadwire on the Market

60" lead wire. Gives the clinician more flexibility while positioning the defibrillator. Is a great benefit with larger patients.

Color Keyed Pouches

Eliminate mistakes—Medtronic is in red, Philips in green, Zoll in yellow, HP in orange and the Anderson connector is blue.

Peel Pouch Bags

Easier for clinicians to open. Pouch also has traditional tear notch.

User Instructions on Release Liner

User instructions and check list on the colored release liner. Colored liner does not become a slipping hazard for the staff on the treatment floor.

Smaller Pad Size

Increasing the conductor plate while reducing the overall size of the electrode. Our aggressive conductive hydrogel reduces the need for a large foam border to anchor the pad to the patient.

Color Keyed Pads

Anterior electrode is blue foam and the posterior electrode is white foam. Helps train the staff to look for colors to aid in pad placement.

Instructions on the Package

Pouch instructions designed only for the US market. Complete user instructions, cautions and warnings on the label, not an IFU in the box.