

MEGAMOVER[®] PORTABLE TRANSPORT UNIT

Our innovative patient rescue and transport units are designed to assist in situations where stretchers cannot.

- Multi-positioned handles for added safety and comfort
- Durability allows units to hold weights ranging from 350-1000 lbs.
- Folds away easily – requiring minimal storage space
- Nonwoven/Poly material acts as a fluid barrier
- Cost effective and disposable

MegaMover[®]

- More durable and easier to grip than a flat sheet
- Each unit weighs only one pound

MegaMover[®] Select™

- Multiple grip options allow for safer patient lifts
- Pull straps for situations where pulling a patient is safer than lifting

MegaMover[®] Transport Chair

- Allows for patients to be carried in the seated position
- Eight handles for ergonomic lifting

MEGA MOVER[®]
portable transport units



The MegaMover[®] Transport Chair is ideal for carrying patients in areas of limited space.



Also available with PowerGrips™ for added comfort and durability.



MegaMover[®] Plus accommodates standard backboards for immobilization. *Backboard not included.*



Product Ordering Information

CODE	PRODUCT NAME	MATERIAL	COLOR	SIZE	PACKING
51926	MegaMover [®]	Nonwoven/Poly	White/Black	40" x 80"	10/Case
53376	MegaMover [®] Plus	Nonwoven/Poly	Gray/Black	40" x 80"	10/Case
56307	MegaMover [®] Transport Chair	Nonwoven/Poly	White/Yellow	36" x 36"	10/Case
52037	MegaMover [®] 350l	Nonwoven/Poly	Gray	50" x 84"	25/Case
65459	MegaMover [®] Plus with PowerGrips™	Nonwoven/Poly	Gray/Black	40" x 80"	8/Case
66563	MegaMover [®] with PowerGrips™	Nonwoven/Poly	White/Black	40" x 80"	8/Case
66564	MegaMover [®] Transport Chair w/ PowerGrips™	Nonwoven/Poly	White/Yellow	36" x 36"	8/Case
72657	MegaMover [®] Select™ Transport Unit	Nonwoven/Poly	White/Green	40" x 80"	8/Case