

STOPS SEVERE BLEEDING

CELOX-A™

Pre-filled Applicator





The effective way to treat penetrating trauma

Penetrating wounds like shrapnel or knife wounds are very difficult to treat. They can often be lethal. It is extremely difficult to get a haemostatic agent to the real source of bleeding. Now there is a solution.

Celox-A™ is a unique applicator delivery system designed to give you an easy, safe and effective way to stop life-threatening bleeding from difficult to treat penetrating trauma. The instinctive applicator allows you to get the Celox™ granules through a small entry wound, directly to the bleeding site in just a few seconds. In clinical testing Celox-A™ has repeatedly shown itself to be

able to quickly and reliably stop deep arterial bleeds. Of course that should not surprise you - Celox™ granules have already been shown to result in 100% survival in independent US marines' severe bleeding tests**. Celox's™ unique, patented and natural formula, works independently of the body's normal clotting mechanisms to robustly clot blood even where normal blood clotting is slow or impaired.

Are You Prepared...?




CELOX is suitable for...

- Small entry wounds
- Penetrating wounds
- Knife and shrapnel wounds
- Bullet entry wounds
- Applying through strong blood flows
- Controlled application to all wounds

CELOX has been shown to:

- Save lives
- Reliably stop arterial bleeding**
- Significantly reduce blood loss**
- Clot blood containing anticoagulants like Coumadin (warfarin)*
- Quickly clot hypothermic (cold) blood*

How does Celox work?



Pouring the CELOX into the wound prevents blood loss by forming a gel-like plug as the CELOX links to the surface of red blood cells.

Product ordering information

Order codes:

Box (70): FG08832021

Box Size: 170mm x 330mm x 390mm



REGULATORY APPROVAL FOR SALE
IN USA & EUROPE

Medtrade Products Limited

Electra House

Crewe Business Park, Crewe,
Cheshire, UK. CW1 6GL

Tel: (44) 1270 500019

Fax: (44) 1270 500045

sales@medtradeproducts.co.uk

* Details of these tests are available at: www.celoxmedical.com

** Journal Academic Emergency Medicine Jan 08, Pages 74-81