# Flexi-Slip<sup>™</sup> Stylet

### RUSCH







## **ET Tube Guide**

The Rüsch Flexi-Slip™ endotracheal tube stylet has a low friction coating to facilitate smooth withdrawal. It is designed with a soft distal tip to reduce the risk of trauma. Excellent malleability allows easy shaping and forming. Also available in a set, pre-loaded inside an endotracheal tube to save time and storage space.

#### Flexi-Slip<sup>™</sup>

- Soft distal tip reduces risk of trauma
- Low friction coating for easy removal
- ♦ Packaged sterile
- ♦ Single use device
- ♦ Latex-free

#### Flexi-Set™

- ♦ Safety Clear Plus™ Murphy eye endotracheal tube
- ♦ Pre-assembled and ready to use
- Available in a full range of sizes cuffed and uncuffed

Fiexi-Siib	- E I	Tube Guia	е
Dun dunat #	C:	Fito Tulos	7

			_	
Product #	Size	Fits Tube	Length	Box Qty
502501	6Fr	2.0-3.5mm	30cm	20
502503	10Fr	4.0-6.5mm	44cm	20
502505	12Fr	5.0-10mm	45cm	20
502507	14Fr	7.0-10mm	45cm	20

Flexi-Set™ - Pre-loaded Uncuffed ET Tube				
Product #	Size/I.D	Box Qty		
506525	2.5mm	10		
506530	3.0mm	10		
506535	3.5mm	10		
506540	4.0mm	10		
506545	4.5mm	10		
506550	5.0mm	10		

Lit No: L1008-0100

Product #	Size/I.D.	Box Qty
504550	5.0mm	10
504555	5.5mm	10
504560	6.0mm	10
504565	6.5mm	10
504570	7.0mm	10
504575	7.5mm	10
504580	8.0mm	10
504585	8.5mm	10
504590	9.0mm	10

To place an order contact
Customer Service at **1-800-514-7234** 

