# Measure to Improve



# **CPR**meter 2

In the early critical moments of a Cardiac Arrest, the CPRmeter 2 helps first responders optimise CPR performance by providing real-time feedback on essential parameters of CPR. Immediately afterwards, the Quick Review event statistics enable the responding team to quickly self-evaluate their performance and provide valuable insights for improvement next time.



## **CPR**meter 2



### Specifications

#### PHYSICAL

Weight

162g (5.7 oz)

Size HxWxD

153 mm x 64 mm x 25 mm

28 mm x 35 mm

128 x 160 pixels

Flashing green light

Solid Amber light

Spanner icon

TFT display

Ruggedness

### IP55 and 1 meter drop test

#### OPERATION

**Display** Dimensions Type Resolution

#### Display indicators

Battery OK Low Battery

Device Error Warning Service Required

#### **CPR** Targets

Compression Depth Compression Release Target

Compression Rate Target Compression Counter Inactivity timer ≥ 50 mm (2") ±10 %
< 2.5 kg (5.5 lbs) force: + 1.5kg to -2.0 kg (+ 3.3 lbs to 4.4 lbs)</li>
100 to 120/min ± 3/min
1-999 (reset after 5 secs)
0.03-9.59 (min:sec)

Small Low battery icon on screen (turning on)

Large Low battery icon on screen (turning off)

## Ordering information

Part No. 801-00249 801-00243 801-00240 Description CPRmeter 2 (Multi-Bundle) CPRmeter 2 (AP Bundle) CPRmeter 2 (Nordic)

Patient Adhesives Battery Hatch Replacement (3 pk)

Accessories 801-10150 801-10550

Consumables 801-10850

801-10360

Red Sleeve CPRmeter 2 Rugged Carry Case

#### Type Capacity Standby life

**Data** Storage Transfer 300 minutes of data or 20 CPR Sessions

2 years (after two years min 30mins CPR)

minimum 10 episodes of 30min CPR

2 x 1.5V AAA

Bluetooth Smart

EMC CLASSIFICATION Meets IEC 60601-1-2 and RTCA/DO-160F

#### ENVIRONMENTAL CONDITIONS

Storage temperature Relative humidity Operating temperature -20° to 70° C (-4° to 158° F) 5% to 75% 0° to 50° C (32° to 122° F)

#### Patient Adhesives

Size Material

Shelf life

39mm x 90mm (1.5" x3.5") Foam pad with biocompatible adhesive on each side. 2 years when applied to CPRmeter or 4 years unopened packaging

### Warranty

One Year limited. Please see Laerdal Global Warranty for terms and conditions

Ergonomic design

Ē

AAA Batteries

 ►
 **↗** 

 ►
 ►

 ►
 ►

 Bigger Screen



Low cost of ownership

+

Easy to clean

A Bluetooth Smart



For more information, visit www.laerdal.com

022017V02 @2017 Laerdal Medical. All rights reserved.