

smiths medical
jelco®




**ViaValve®
Winged Safety
IV Catheter with
Blood Control**

Designed to protect against
blood exposure, contamination
and needlestick injury



VIAVALVE® WINGED SAFETY IV CATHETERS

TRUST IS EVERYTHING



ORDERING INFORMATION

ViaValve® Winged Safety IV Catheters					
ViaValve® Winged Codes	Gauge Size / Length /Catheter Material	Features	Flow Rate mL/min	Units Per Case	
3282	16G x 1 1/4" Tehrilon® Polyurethane	Winged Hub, Radiopaque	220	200	
3285	18G x 1 1/4" Tehrilon® Polyurethane	Winged Hub, Radiopaque	110	200	
3287*	20G x 1" Tehrilon® Polyurethane	Winged Hub, Radiopaque	65	200	
3286*	20G x 1 1/4" Tehrilon® Polyurethane	Winged Hub, Radiopaque	63	200	
3280*	22G x 1" Tehrilon® Polyurethane	Winged Hub, Radiopaque	38	200	
3283*	24G x 5/8" Tehrilon® Polyurethane	Winged Hub, Radiopaque	24	200	

*Features FLASH-VUE® notched needle for early flash detection

References

- Differences in penetration force of intravenous catheters: effect of grinding methods on inner needles of intravenous catheters. Suzuki T, Tanaka A, Fukuyama H, Nishiyama J, Kanazawa M, Oda M, Takahashi M. Tokai J Exp Clin Med. 2004 Dec;29(4):175-81.

To learn more about ViaValve® safety IV catheters and schedule a demonstration, contact your local Smiths Medical sales associate or visit www.smiths-medical.com

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES

Smiths Medical ASD, Inc.
 6000 Nathan Lane North
 Minneapolis, MN 55442, USA
 Tel: 1-614-210-7300
 Toll-Free USA: 1-800-258-5361
www.smiths-medical.com

MSOUC-1069

**Rx
 ONLY**

smiths medical

Find your local contact information at: www.smiths-medical.com/customer-support

Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. ViaValve, ViaValve Winged, FLASH-VUE, Tehrilon and the Smiths Medical and Jelco design mark are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trademarks or service marks of their respective owners. ©2016 Smiths Medical. All rights reserved. SS195431EN-082016