



### 12 LEAD ECG KIT

#### V1 - V2:


Small rectangular electrode with two sensors

#### V3 - V6:

Small rectangular shaped electrode for placement in tight areas

#### RA, LA, RR, LL:

Large round electrodes for superior adhesion



### SINGLE USE PACKS

Packaged in no-waste quantities



Quantity clearly defined on pouch for ease of use in emergency situations




# EMS Electrodes


## 12 Lead ECG Kit

Built with time savings in mind, our single use pouch contains 10 electrode perfect for EMS diagnostic stress tests.

### A10052-10S - STRESS TEST KIT - FOR 10 LEAD PLACEMENT 10 ELECTRODES/ KIT | 50 KITS / CASE

Size	Shape	Per Pouch	Wet Gel
2.25" / 57.15mm	Round	4	
1.875" x 1.25" / 47.625mm x 31.75mm	Rectangle	4	
1.75" x 1.625" / 44.43mm x 41.275mm	Square, dual element	1	

## WHAT'S IN A KIT?



**4** 2.25" / 57.15mm  
Foam | Wet Gel

**4** 1.875" x 1.25" /  
47.625mm x 31.75mm  
Foam | Wet Gel

**1** 1.75" x 1.625" /  
44.45mm x 41.275mm  
Perforated | Foam | Wet Gel



Perforation allows separation when desired