

# InfuseIT™

DISPOSABLE PRESSURE INFUSOR

500 ml

1000 ml

3000 ml



CE 0086

**Zefon**  
INTERNATIONAL

5350 SW 1st Lane  
Ocala, FL 34474  
Tel: +1352-854-8080  
Fax: +1352-854-7480  
www.zefon.com

EC REP Authorized European Representative

Blue Box Medical Limited  
Unit 29, New Forest Enterprise Centre,  
Rushington Business Park, Chapel Lane  
Totton, Southampton, SO40 9LA  
Tel: 023 8066 9000; Fax: 023 8066 9002

## Applications:

- Pre-Operative Care
- Labor & Delivery
- Radiology; Special Procedures
- Emergency/Trauma Department
- Cardiac Catheterization Laboratory
- Post-Operative Care
- Intensive Care
- One Day Surgery
- Surgery

### Features:

- Unique pressure gauge is easy to read and durable.
- Convenient IV pole loop hanger is sturdy and easy for personnel to hang onto IV poles regardless of height.
- IV fluid bag hook is durable and designed to fit all standard bags snugly.
- Durable Construction - will not rip or tear.
- Clear sleeve around bag provides a visible window to see what fluids are being pressurized.
- Palm fitted, efficient hand bulb makes inflation easy.
- Choice of stopcock valve or bulb thumbwheel valve.



### Easy to Read Pressure Gauge

Unique gauge provides an easy to read 360 degree window, which can be read from far away.



#### Starting Position

Displays the pressure at 0mmHg before any pressure is applied.



#### Normal Pressure

Indicator marks are located at the 175mmHg & 300mmHg levels. Green color at 300mmHg indicates normal proper inflation.



#### Over Pressurization

Red color indicates over pressurization. When pressure reaches this area an automatic bleed valve will engage reducing the pressure to a safe level.

Product Number	Description
ZIT-500	InfuseIT 500ml with Stopcock valve, 10/box
ZIT-502	InfuseIT 500ml with Thumbwheel valve, 10/box
ZIT-1000	InfuseIT 1000ml with Stopcock valve, 10/box
ZIT-1002	InfuseIT 1000ml with Thumbwheel valve, 10/box
ZIT-3000	InfuseIT 3000ml with Stopcock valve, 10/box
ZIT-3002	InfuseIT 3000ml with Thumbwheel valve, 10/box

Distributed by: