

Microflex®

FreeForm® SE

Powder-Free Nitrile Exam Gloves

Item

Extra Small.....	FFS-700-XS
Small	FFS-700-S
Medium.....	FFS-700-M
Large.....	FFS-700-L
Extra Large.....	FFS-700-XL



Recommended for use in the following industries:

Dental / EMS / Healthcare / Industrial Safety / Laboratory

If your need is softness combined with durability, **FreeForm® SE** is the glove for you. Designed with strength and comfort in mind, **FreeForm® SE** gloves are a unique combination of features in a nitrile glove.

FreeForm® SE gloves are the best overall value in a comfortable mid-range synthetic glove.

Recommended for use when:

- Working in wet environments that demand a reliable grip
- Light to mid-duty tasks require a combination of comfort and durability
- Precise contact and softness require enhanced tactile sensitivity

Microflex® FreeForm® SE Profile

For current specification information, please visit www.microflex.com.

	ASTM D 6319	FreeForm® SE		
Length (mm/in):	230 / 9.1	243 / 9.6	Material: Nitrile	Type: Non-sterile
Cuff Thickness (mils):	n/a	2.8	Size: XS - XL	Color: Blue
Palm Thickness (mils):	min. 2.0	3.5	Shape: Ambidextrous	Usage: For single use only
Finger Thickness (mils):	min. 2.0	5.1	Cuff: Standard examination length (9.5"), Beaded	
Tensile Strength:			Glove Exterior: Textured Fingers	
Before Aging (MPa):	min. 14	32	Glove Interior: Powder-Free, Chlorinated	
After Aging (MPa):	min. 14	31	Packaging: 1,000 gloves/case, by weight; 10 boxes/case	
Elasticity/Elongation:			100 gloves per dispenser box, by weight	
Before Aging (%):	min. 500	500	Caution: Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.	
After Aging (%):	min. 400	400		

The reason is the quality... the difference is Microflex®.

See for yourself why Microflex® gloves are the first choice of exam glove users nationwide. For **FREE SAMPLES**, please visit www.microflex.com. For the name of a distributor in your area, call: **(800) 876-6866**.

MICROFLEX®
THE MOST TRUSTED NAME IN GLOVES®
A division of BarrierSafe Solutions International, Inc.