



Put the Science of Protection to Work in Your Lab

KIMBERLY-CLARK* PURPLE NITRILE* Exam Gloves

Your employees are your most valuable asset. That's why they deserve the worry-free protection of KIMBERLY-CLARK* PURPLE NITRILE* Exam Gloves. Our gloves guard against incidental splash exposure. The gloves are made of a highly durable nitrile polymer that eliminates the risks and costs associated with type 1 natural rubber latex allergies.



KIMBERLY-CLARK* PURPLE NITRILE* Exam Gloves are the gloves of choice in leading labs because they are comfortable to wear and provide excellent tactile sensitivity. When users

see the distinctive color purple, they know they can concentrate on the task at hand, protected from known and unknown risks.



KIMTECH*
BRAND



KIMBERLY-CLARK® PURPLE NITRILE® Exam Glove Features and Specifications

- Latex-free
- Powder-free
- Beaded cuff
- 6 mil thickness
- Ambidextrous
- Textured fingertips
- 9.5" & 12" lengths
- AQL 1.5
- Excellent dexterity and tactile sensitivity
- Comfortable fit
- KIMBERLY-CLARK® PURPLE NITRILE-XTRA® material was tested for use with 10 chemotherapy drugs and glutaraldehyde

Product Specifications

Gauge Thickness Measurements	MM	MIL
Middle Finger	.15	5.9
Palm	.12	4.7
Cuff	.09	3.5
Average Length	242mm (9.5") XTRA® 305mm (12.0")	

Physical Properties

Before Aging	
Tensile Strength	21 MPa
Ultimate Elongation	550%
After Aging	
Tensile Strength	21 MPa
Ultimate Elongation	500%

Quality Standards

Exceeds current ASTM 6319 standard for critical defects (AQL-2.5). AQL for critical defect is 1.5.
Manufactured in accordance with Quality System ISO 9001:2000.

Best Protection

When you're in the market for protection, turn to the proven leader.
KIMBERLY-CLARK® PURPLE NITRILE® Exam Gloves provide peace of mind for your lab personnel.



KIMBERLY-CLARK® PURPLE NITRILE® Exam Gloves

Description	Color	X-Small	Small	Medium	Large	X-Large	Gloves/Boxes	Total/Case
KIMBERLY-CLARK® PURPLE NITRILE® Exam Gloves 9.5" Ambi	Purple	55080	55081	55082	55083	55084	100/10 (90/10)	1000 (XL = 900)
KIMBERLY-CLARK® PURPLE NITRILE-XTRA® Exam Gloves 12" Ambi	Purple	55090	50601	50602	50603	50604	50/10	500

Specific chemical resistance data can be found at: www.kc-safety.com/chemicalbarrierdata

Our Guarantee

Your total satisfaction means everything to us. If, for any reason, our products do not meet your expectations, Kimberly-Clark will reimburse you† for your initial purchase, via FREE product, for up to \$1,000.

For more information on KIMBERLY-CLARK PROFESSIONAL®, visit us online at www.kcprofessional.com, contact your Kimberly-Clark Sales Representative, or call us at 1-888-346-GOKC (4652).

† Guarantee extended to consuming end-user accounts only.

Kimberly-Clark warrants that its products comply with K-C's standard specifications as of the delivery date to K-C's authorized distributors/direct purchasers. Except to the extent prohibited by applicable law: 1) **THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE;** 2) K-C is not liable for any kind of special, incidental, or consequential damages; and 3) K-C's liability for breach of contract, tort, or other cause of action shall not exceed the product purchase price. Purchasers and users are deemed to have accepted the above warranty and limitation of liability, and cannot change the terms by verbal agreement or by any writing not signed by K-C. To the extent required by applicable law, K-C does not limit its liability for death/injury resulting from K-C's negligence.

®/* Trademarks of Kimberly-Clark Worldwide, Inc. Marques déposées de Kimberly-Clark Worldwide, Inc. The COLOR PURPLE is a Registered Trademark of KCWW. © 2009 KCWW.
Printed in the USA. K01716 K3512-09-01 10/09

Effective October 2009