

A.C.E  
Airway Circulation Evaluator

# STAT-Check

## CO2 Indicator & Resuscitation Bag

**Latex Free**



**IT'S SAFE TO BE SURE !!!**

The A.C.E. **STAT-Check** indicator is not located on the exhalation port, therefore reducing risk of ambient CO2 indicating a false positive.

- Integrated Colorimetric CO2 Indicator
- Small Compact Design
- Simple to use: Easy to view, one color Interpretation
- A.C.E. **STAT-Check** also sold separately

Stand-Alone A.C.E. **STAT-Check** Indicator



ACE/ **STAT-Check** Indicator Arm  
Exhalation Port



Patient Connection Port      Double swivel Elbow

Parameters	Adult	Pediatric
<b>Bag Volume</b>	<b>1900ml</b>	<b>1000ml</b>
<b>Stroke Volume (ASTM Standard Hand)</b>	<b>&gt;960ml</b>	<b>450-700 ml</b>
<b>F<sub>D</sub>O<sub>2</sub></b>	<b>99%</b>	<b>99%</b>
<b>BPM (@ room temp)</b>	<b>&gt;65</b>	<b>&gt;65</b>
<b>Reorder Numbers</b>	<b>SC9000 series</b>	<b>SC8000 series</b>

Manufactured by:  
**Ventlab Corporation**  
155 Boyce Drive  
Mocksville, N.C. 27028  
Toll-Free: 800-593-4654  
Phone: 336-753-5000  
Fax: 336-753-5002

Infant Resuscitation bags available with  
“stand-alone” **STAT-Check**

US Patent 4945918 and others pending