



Endotracheal Tube Guide

The Rüscher Slick® Stylet is made of malleable coated wire that holds its shape to facilitate intubation. The specially formulated nylon coating does not crack or peel while allowing for smooth tube insertion. The stylet is fitted with a stopper that helps maintain the desired stylet position within the endotracheal tube. For added convenience and economy of space, Slick® Stylet is available pre-assembled in an endotracheal tube.

Slick® Stylet

- ✧ Available in 3 color coded sizes for easy identification
- ✧ Handle forming capability for a better grip
- ✧ Sterile, single use device
- ✧ Malleable coated wire
- ✧ Latex-free

Slick Set®

- ✧ **NEW!** Corkstyle topper helps maintain stylet position
- ✧ Pre-assembled and ready to use - Saves you both time and space
- ✧ Latex-Free
- ✧ Available in a full range of sizes - cuffed and uncuffed
- ✧ Safety Clear Plus™ Murphy eye endotracheal tube

Slick® Stylet - ET Tube Guide

Product #	Size	Color	Fits Tube	Box Qty
500	6Fr	Blue	2.0 mm - 3.5 mm	25
750	8Fr	Grey	4.0 mm - 6.5 mm	25
1000	10Fr	Green	7.0 mm - 10.0 mm	25

Slick Set® - Pre-loaded Cuffed ET Tube

Product #	Size/I.D.	Box Qty
170050	5.0 mm	10
170055	5.5 mm	10
170060	6.0 mm	10
170065	6.5 mm	10
170070	7.0 mm	10
170075	7.5 mm	10
170080	8.0 mm	10
170085	8.5 mm	10
170090	9.0 mm	10
170095	9.5 mm	10
1700100	10.0 mm	10

Slick Set® - Pre-loaded Uncuffed ET Tube

Product #	Size/I.D.	Box Qty
150020	2.0 mm	10
150025	2.5 mm	10
150030	3.0 mm	10
150035	3.5 mm	10
150040	4.0 mm	10
150045	4.5 mm	10
150050	5.0 mm	10
150055	5.5 mm	10

Latex-Free

To place an order contact
Customer Service at **1-800-514-7234**