



Bariatric Transfer Flat

Introducing the Bariatric Patient Transfer Flat

There are many patient transfer flats available on the market today. None of them qualify as a spinal immobilization device. Unfortunately, the patient who typically requires this device is found on the ground from a fall and requires spinal immobilization.

The patented technology of the Thomas Transfer Flat is designed to accommodate most backboards and qualifies as a spinal immobilization device.

NO MORE LIFTING THE PATIENT TWICE!



Your backboard can simply be attached to the center of the Thomas Transfer Flat with access to pins and hand holds for strapping. The ten outside ergonomic handles allow for even distribution of carrying.

The Thomas Transfer Flat is designed to roll up and be self contained. Simply unbuckle, unroll, slide in the backboard* and your ready to transport. The heavy-duty vinyl construction is easily cleaned with soap and water - no Velcro. Designed by EMTs, Paramedics and Physicians, the Thomas Transfer Flat is ready for your real needs in the field.

Specifications:

Length: 72" Width: 46" Weight: 7 lbs
 Maximum Weight Capacity: 1000 lbs
 Operators: 2 -10

*Backboard not included

Item #	Description	Price
TTS200	Bariatric Transfer Flat	\$178.50

Thomas EMS
 PO Box 651305
 Salt Lake City, Utah 84165
 Toll Free (800) 445-3640
 (801) 262-6503 • Fax (801) 268-9272
www.thomasems.com