

Microflex®

ComfortGrip®

Powder-Free Latex Exam Gloves

Item

Extra Small..... CFG-900-XS
Small CFG-900-S
Medium..... CFG-900-M
Large..... CFG-900-L
Extra Large..... CFG-900-XL



Recommended for use in the following industries:

Automotive / Dental / EMS / Healthcare / Industrial Safety / Laboratory

If you're looking for the best combination in value and quality, **ComfortGrip®** is the latex glove you need. **ComfortGrip®** gloves are backed by the industry's leading quality system to offer consistent comfort and quality in every box. Fully textured for a secure grip, **ComfortGrip®** gloves provide the fit and protection you need at a competitive price.

ComfortGrip® gloves offer the best overall value and quality in a latex glove.

Recommended for use when:

- A secure grip is needed in wet and dry environments
- Working with small instruments requires maximum sensitivity
- Comfort, quality, and value are crucial to your needs

Microflex® ComfortGrip® Profile			For current specification information, please visit www.microflex.com .	
	ASTM D 3578	ComfortGrip®		
Length (mm/in):	230 / 9.1	245 / 9.6	Material: Latex	Type: Non-sterile
Cuff Thickness (mils):	n/a	3.9	Size: XS - XL	Color: Natural
Palm Thickness (mils):	min. 3.1	5.1	Shape: Ambidextrous	Usage: For single use only
Finger Thickness (mils):	min. 3.1	5.9	Cuff: Standard examination length (9.5"), Beaded	
Tensile Strength:			Glove Exterior: Fully Textured, Polymer-Coated	
Before Aging (MPa):	min. 18	23	Glove Interior: Powder-Free, Chlorinated	
After Aging (MPa):	min. 14	19	Packaging: 1,000 gloves/case, by weight; 10 boxes/case	
Elasticity/Elongation:			100 gloves per dispenser box, by weight	
Before Aging (%):	min. 650	650	Caution: This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.	
After Aging (%):	min. 500	500		

The reason is the quality... the difference is Microflex®.

See for yourself why Microflex® gloves are the first choice of exam glove users nationwide. For **FREE SAMPLES**, please visit www.microflex.com. For the name of a distributor in your area, call: **(800) 876-6866**.

MICROFLEX®
THE MOST TRUSTED NAME IN GLOVES®
A division of BarrierSafe Solutions International, Inc.