


# QD-4 Traction Leg Splint



**Comes in complete kit with straps, ankle hitch and carrying case included.**

**Hex Lock allows for easy adjustment, in the worst of conditions.**



**Hard Anodized Aluminum Ratchet**



**Injection Molded Ratchet**


- ⊕ Compact, Strong, Light
- ⊕ QD-3: Pediatric Unit
- ⊕ Aluminum or Molded Ratchet
- ⊕ Integral Strap System
- ⊕ Field tested
- ⊕ QD-4 unit fits PEDI to Adult
- ⊕ Minimal patient movement
- ⊕ Intuitive Adjusting System

Centering the Hex Lock allows for quick length adjustment. (1) Turning the Hex Lock in either direction locks the splint to the desired size. (2)

1.)



2.)



## **Labeled Carrying Case**



**FareTec Inc.**  
Fracture and Rehabilitation Equipment

1610 West Jackson St #6  
Painesville OH 44077

440 350 9510 v  
440 350 9520 fax  
800 322 3273

www.faretec.com  
Service@faretec.com

