

THOMAS ALS HARD CASE

Innovative Solutions With You In Mind



Keep your ALS equipment and pharmaceuticals in a lockable, durable and organized hard case.

Large Removable Drug Case

The Thomas ALS Medication Hard Case features a Large Removable Padded Drug Case located in the lid of the Case. This allows several medics the flexibility to work out of the same Case in multiple victim situations. In single victim situations, the Drug Case can remain secured to the lid and still fully functional.

The Padded Drug Box is secured to the lid of the Hard Case by 6 neodymium magnets, completely eliminating the need for Velcro. The use of fluid-resistant heavy weight vinyl and closed cell foam also significantly reduces breakage risk and makes decontamination and basic cleaning quick and easy.

The layout of the elastic drug securing system and pocket design allows for the accommodation of all methods of drug delivery and storage (pre-filled, ampoules, boxed storage, etc..) The Drug Case also contains a clear vinyl locking section for narcotic storage.



Foam Divided Base

The base of the Thomas ALS Medication Hard Case consists of a modular, foam divider system plus 2 removable fold outs featuring clear vinyl, color coded pockets. The adjustable foam dividers have been shortened in height slightly, to lessen the depth of the base section and to accommodate the two fold out sections.

The base of the hard case can easily accommodate up to 3 liters of IV solution and accessories, your IO device, as well as additional medications. This design allows for all base contents to be quickly and easily viewed upon the opening of the Case. No more digging through the base of your Hard Case searching for equipment.

Interior Dimensions: 18.43" x 14" x 7.62"



Neodymium magnets maintain 99.5% of their strength for up to 100 years!!

Item	Description	Price
TT424	ALS Hard Case Complete, Orange	\$353.60
TT423	ALS Hard Case Complete, Gray	\$353.60
TT434	ALS Hard Case Complete, Yellow	\$353.60
TT321	Large Removable Padded Drug Case	\$135.00
TT742	Modular Foam Divider System	\$60.00

Complete ALS Hard Case includes TT321, TT742 and the Pelican 1550 Case.



THOMAS AIRWAY HARD CASE

Large Intubation Kit

The Thomas Airway Hard Case features a large removable Intubation Kit, located in the lid of the Case. The Intubation Kit can be removed and used by several medics in a multiple victim scenario, or remain secured to the lid and still fully functional. The Intubation Kit is secured to the lid by 6 neodymium magnets, which completely eliminates the need for Velcro on the outside of the kit.

The Intubation Kit is designed to hold multiple intubation handles, Macintosh and Miller blades, ET Tubes, Magill forceps, airways, RSI medications, and all other needed intubation accessories. The use of fluid-resistant heavy weight vinyl on the outside of the intubation kit makes decontamination easy.

Foam Divided Base

The base of the Thomas Airway Case consists of a modular, foam divider system plus 2 removable fold outs featuring clear vinyl, color coded pockets. The base of the Airway Case is designed to accommodate multiple adult and pediatric BVM's.

Item	Description	Price
TT431	Airway Hard Case Complete, Orange	\$357.50
TT432	Airway Hard Case Complete, Gray	\$357.50
TT433	Airway Hard Case Complete, Yellow	\$357.50
TT321	Large Removable Padded Drug Case	\$135.00
TT742	Modular Foam Divider System	\$60.00

Complete Airway Hard Case includes TT341, TT742 and the Pelican 1550 Case.



Thomas BLS Hard Case BLS Lid Section

The Thomas BLS Hard Case features a removable lid section that contains 6 color coded pockets, designed to hold all of your BLS response equipment. The BLS features the same magnetic lid securing system as the other Thomas Hard Cases, and is made out of the same fluid-resistant material.

Foam Divided Base

The base of the Thomas BLS Hard Case consists of a modular, foam divider system plus 2 removable fold outs featuring clear vinyl, color coded pockets. This design allows for all base contents to be protected and quickly and easily viewed upon the opening of the Case. The base of the Thomas BLS Hard Case can hold all of the contents that the lid can accommodate, thus allowing for two EMT's to work on separate patients simultaneously from the same BLS Hard Case if needed.

Item	Description	Price
TT428	BLS Hard Case Complete, Gray	\$345.00
TT429	BLS Hard Case Complete, Orange	\$345.00
TT430	BLS Hard Case Complete, Yellow	\$345.00
TT321	Large Removable Padded Drug Case	\$128.50
TT742	Modular Foam Divider System	\$60.00

Complete BLS Hard Case includes TT331, TT742 and the Pelican 1550 Case.

THOMAS BLS HARD CASE



ALL THOMAS HARD CASE INSERTS ARE DESIGNED FOR USE IN THE PELICAN 1550 OR EQUIVALENT SIZED CASE.