



# AED Pack



## with External Removable Padded Pocket



The Thomas AED pack is designed to safely secure, transport, and store virtually all models of automatic external defibrillators on the market today.

The outside of the pack is clearly marked for quick and easy identification of the unit and the large front external pocket is capable of carrying even the largest AED currently in use. The AED pocket is lined with 1/2 inch Evazote foam on all sides with to keep the AED safe during any usage.

**Inside the pack**, a redi-splint on all sides, which can also be removed and used in any situation that requires splinting. The large clear vinyl pocket located immediately behind the AED allows for the storage of multiple adult electrodes. The location and size of the vinyl pocket also allows the electrodes to be stored in the "attached" position, saving valuable time at the scene.

**Three additional vinyl pockets** are also provided for the storage of pediatric electrodes, clippers, gloves, or razors. There is more than enough room for a myriad of other supplies, and general BLS gear.

**Manufactured with 1000 Dupont denier coated Cordura** and 400D pack cloth, YKK reverse pulls and zipper, The Thomas AED pack also contains multiple D rings for quick and easy attachment of backpack straps or attachment to other packs.



Item	Description	Price
TT831	AED Pack with External Removable Padded Pocket	\$170.00

**Thomas EMS**  
Toll Free (800) 445-3640  
(801) 262-6503 • Fax (801) 268-9272  
[www.thomasems.com](http://www.thomasems.com)