

# HYFIN® CHEST SEAL *Individual Occlusive*

REV101917



## Specifications:

- Device Dimensions:  
6 in. x 6 in.  
(excluding tab)
- Pkg Dimensions:  
H 4.5 in. x W 7.75 in  
x D 0.1 in.
- Weight: 1.31 oz.
- U.S. Patents:  
Patent(s) Pending
- Includes No  
Vent Channels

The HyFin® Chest Seal is a self-adhesive occlusive dressing designed for the treatment and management of an open chest wound caused by a penetrating trauma, such as a gunshot or stab wound. A penetrating chest wound creates an opening in the chest wall through which air can enter the chest cavity; commonly referred to as an “open pneumothorax.” The HyFin® Chest Seal utilizes an enhanced adhesive gel formulation designed to create a superior seal even in extreme conditions to treat and manage an open pneumothorax while preventing the development of a “tension pneumothorax,” which is one of the top three leading causes of preventable deaths on the battlefield. The HyFin® Chest Seal is specifically designed to overcome adverse conditions in austere environments where the casualty may be covered in blood, sweat, body hair, or other environmental contaminants when a reliable chest seal is mandatory.

Featuring a large 6 x 6 in. size, the transparent design allows for accurate placement directly over the wound, easily conforming to the casualty’s chest. A large, easily-visible red elliptical tab facilitating single step peel-and-adhere application is incorporated into the HyFin®’s design. Made in America and 100% latex free, each chest seal is packaged using its signature Red Tip Technology® red tear notches in a rugged, easy-to-open, foil pouch to ensure high quality performance.

## Features:

- Completely occlusive dressing designed for the treatment and management of an open chest wound caused by a penetrating trauma
- New, more compact packaging
- Advanced adhesive technology designed to create a superior seal
- Transparent backing designed to allow placement directly over the wound and easily conform to the casualty’s chest
- Red elliptical tab for single step peel-and-adhere application that can also facilitate burping of the wound

### HyFin® Individual Occlusive Chest Seal

Item #	NSN#
10-0015	6510-01-532-8019



**NORTH AMERICAN RESCUE®**  
www.NARescue.com • 888.689.6277

**PRODUCTS**  
WITH  **MISSION®**