



The perfect fit for neutral alignment

The Ambu Perfit collars are designed to assist with the maintenance of neutral alignment, prevention of lateral sway and anterior-posterior flexion and extension of the spine during transport and routine patient care or movement.

Ambu® Perfit® Extrication Collars

All Ambu Extrication Collars have the following features

- Ambu collars are fully disposable. You no longer have to worry about collars being returned to the source of origin and cross contamination issues are avoided
- Ambu collars come with a pre-shaped flip-chin piece to ensure they can be stored completely flat. Furthermore the flip-chin piece facilitates intubation and removal of vomit from a patient's mouth
- One-piece extrication collar simplifies application and at the same time there are no extra pieces that can get lost
- Nasal cannula holders are located on the collar body to secure the oxygen tubes (except for infant)
- All collars are radiotranslucent, CAT & MRI compatible and come with a ventilated posterior shell for fluid drainage and a large hole for easy access. All of these features allow for patient examinations to take place safely without removing the collar.

INMOVILIZACIÓN

**AMBU® PERFIT EXTRICATION COLLAR
IS A ONE-SIZE COLLAR**

It comes in four adult sizes, a pediatric size and an infant size.

TECHNICAL DATA

Shell material:	polyethylene
Sizes:	1 (infant), 2 (paediatric), 3 (neckless), 4 (short), 5 (regular) and 6 (tall)
Colour:	White/light grey

ORDERING INFORMATION

	CAT. NO.
Perfit Infant: size 1 / cartons of 20 pcs.	264 201 000
Perfit Paediatric: size 2 / cartons of 20 pcs.	264 202 000
Perfit Neckless: size 3 / cartons of 20 pcs.	264 203 000
Perfit Short: size 4 / cartons of 20 pcs.	264 204 000
Perfit Regular: size 5 / cartons of 20 pcs.	264 205 000
Perfit Tall: size 6 / cartons of 20 pcs.	264 206 000
Carry bag for collars, fits 6 pieces / cartons of 20 pcs.	264 008 000
Sample pack, one of each Perfit size	264 020 000



Ambu A/S

Baltorpbakken 13
DK-2750 Ballerup
Denmark
T +45 72 25 20 00
F +45 72 25 20 50
www.ambu.com

INMOVILIZACIÓN

