



FILTER TECHNOLOGY



# GVS ELIPSE P100 UNVALVED

COMPACT AND LIGHTWEIGHT HALF MASK RESPIRATOR TO PROTECT AGAINST MICROORGANISMS AND VIRUSES, THE UNVALVED VERSION PREVENTS CROSS CONTAMINATION FROM THE WEARER

## FEATURES AND BENEFITS

- ◆ FILTER EFFICIENCY EXCEEDS 99.97% FOR PARTICULATE > 0.3 MICRONS
- ◆ NO VALVE TO PROTECT WEARER AND PATIENT
- ◆ MEDICAL GRADE HYPO-ALLERGENIC THERMO PLASTIC ELASTOMERE
- ◆ LATEX, SILICONE AND ODOR FREE
- ◆ AVAILABLE IN 2 SIZES (S/M, M/L)
- ◆ FILTERS ARE REUSABLE AND REPLACEABLE
- ◆ POSSIBILITY TO DISINFECT MASK AND FILTERS
- ◆ INDIVIDUALLY PACKAGED
- ◆ ADJUSTABLE STRAPS WITH A 4-POINT HARNESS
- ◆ EASY COMMUNICATION
- ◆ MADE IN FINDLAY, OHIO
- ◆ NIOSH APPROVAL PENDING (APRIL 2021)



FILTER TECHNOLOGY

**SPR643**

S/M Elipse P100 Unvalved  
Mask & Filters

**SPR644**

M/L Elipse P100 Unvalved  
Mask & Filters



Replacement P100 Filter  
(SPR321)

Hard carry case  
(SPM001)

TSl portacount fit test adapter  
(SPM414)



**1** CLOSED VALVE



**2** LARGE PLEATED OPEN AREA  
WITH NO VALVE = INSPIRE/ EXPIRE  
THROUGH FILTER