


MERET

MERET

877.222.0200
www.MeretUSA.com
info@meretusa.com
Ph: 714.259.4700
Fx: 714.259.4749




MERET OXYGEN REGULATORS

MERET oxygen regulators incorporate the latest in technology to deliver the highest level of safety, reliability, durability, and precision flow. EMS series unibody regulators are constructed from high strength solid brass with brass internal high pressure flow components to ensure added strength and safe operation.



MEDICAL * EMERGENCY * EQUIPMENT & TECHNOLOGY

MAGNUM SERIES ALL-BRASS OXYGEN REGULATOR



Product Dimensions:
6.5" x 3.5" x 3.5"

Product Weight:
2 lbs 1 oz

MERET CGA REGULATOR OPTIONS

REGULATOR	MODEL	FLOWS (L/MIN)	OUTLETS	CONSTRUCTION	WARRANTY
MAGNUM	EMSMREG8725-xxx	0-25	2 PTO + Barb	Solid 360 Brass, Chrome Plated	LIFETIME
BRUUT	EMSREG8725-xxx	0-25	0-2 PTO + Barb	Solid 360 Brass, Chrome Plated	5 yr
ECONOMY	AREG8725-B2D-BU	0-25	2 PTO + Barb	Aluminum/ Brass Internal, Anodized	5 yr
MINI	AREG87xxME-BU	0-8, 0-15	Barb Only	Aluminum/ Brass Internal, Anodized	5 yr
FLOWMETER	MFMC-25	0-25	N/A	Aluminum/ Brass Internal, Anodized	5 yr
PRESSURE REDUCER	PRREG54-2D	50 PSI Preset	2 PTO	Solid 360 Brass, Chrome Plated	5 yr

BRUUT




ECONOMY




MINI



FLOWMETER



PRESSURE REDUCER



All regulators are individually calibrated and tested 100% for flow and leak. ASTM G175, ISO 13085, ISO 9001:2000