


Decompression Needle



SPECIFICATIONS

- Packaged Size: 5.5" Length x .5" Width
- Needle Size: 14 and 10 Gauge x 3.25 inch
- Unit Weight: ~0.11g
- Twist off cap
- EO Sterilization
- Single Use Only
- Not made with natural rubber latex
- Not made with DEHP

Item #	Description	UOM
MS-81114	Decompression Needle 14G x 3.25"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-87010	Decompression Needle 10G x 3.25"	50 EA/BX, 4 BX/CS, 200 EA/CS