



# MEGA MOVER®

## portable transport chair



The MegaMover® Portable Transport Chair is a strong, disposable alternative for transporting patients down stairs and through narrow hallways or extracting them from confined spaces like bathrooms.

- **Sanitary and disposable.** If it gets soiled, just throw it away.
- Made of fluid-resistant material strengthened with reinforced nylon straps.
- **Strong.** Carries up to 350 pounds with the weight distributed across eight nylon-reinforced handles.
- Carrying handles give a firm grip for up to four people to share the load.
- **Compact.** Weighs only 10 ounces and rolls to the size of a water bottle.



Code	Item	Material	Color	Size	Case	Max Weight
56307	MegaMover® Transport Chair	Nonwoven/Poly	White/Yellow	36" x 36"	1/BG; 10BG/CS	350 lbs.