

REORDER NUMBER: 10244



BARIATRIC TRANSPORT WHEELCHAIR

22" x 18" SEAT

- Designed with sturdy steel frame and dual cross brace for reliable support
- Features foldable back for convenient storage and transport
- Tool-free, swing-away footrests allow for added versatility
- Detachable desk-length armrest
- Equipped with 12-inch rear wheels with push-lock wheel brakes
- Powder-coated finish provides long-lasting protection
- Black-padded nylon upholstery with chart pocket for added convenience
- Limited Lifetime Warranty
- HCPCS: E1039



WEIGHT CAPACITY

- 450 lb.