

# Pulmodyne®

## Control-Cric™

### Control Your Cricothyrotomy.

The Control-Cric System includes the Cric-Knife™ and Cric-Key™. The Cric-Knife has a 10mm long dual sided blade with an integrated sliding tracheal-hook\*. The Cric-Key has a pre-shaped stylet that provides airway confirmation with tactile feel of the tracheal rings, along with a 5.5mm soft tube with cuff.



\*Sliding Tracheal Hook, slides over the dual sided blade and into the incision site, to help you maintain your airway.

Part No.	Product Description	Case/Qty.
313-4564	Control-Cric System (Sterile)	5
313-4565	Control-Cric Training System (Non-Sterile), <b>Not For Human Use</b>	10
313-4567	Control-Cric International Military System (Sterile)	5
313-4568	Control-Cric International Military Training System (Non-Sterile), <b>Not For Human Use</b>	10

For more information about our complete line of ventilation products visit our website at [www.pulmodyne.com](http://www.pulmodyne.com) or call 317.246.5505.

.... bringing change to life®

