



Put the Science of Performance to Work in Your Lab

KIMBERLY-CLARK* STERLING* Nitrile Exam Gloves

KIMBERLY-CLARK* STERLING* Nitrile Exam Gloves bring together the protection of nitrile with the sensitivity of latex and superior strength. That's why they're the top-selling nitrile exam glove in the United States.

Designed for a demanding lab environment, KIMBERLY-CLARK* STERLING* Nitrile Exam Gloves also help reduce environmental impact. Compared to the competitive latex product, you can save 1 ton of waste per year based

on a monthly usage of 25 cases. If you're standardizing to a single exam glove in your lab, the glove of choice is the KIMBERLY-CLARK* STERLING* Nitrile Exam Gloves.

KIMBERLY-CLARK* Green Project
www.kimtech.com/greenproject



STERLING* Nitrile Exam
#1 Glove
in the U.S.



KIMTECH*

BRAND



KIMBERLY-CLARK® STERLING® Nitrile Exam Glove Features and Specifications

- Latex-free
- Powder-free
- Beaded cuff
- 3.5 mil thickness
- Ambidextrous
- Textured fingertips
- 9.5" and 12" lengths
- AQL 1.5
- Provides the protection of nitrile with the sensitivity of latex
- KIMBERLY-CLARK® STERLING® NITRILE-XTRA® was tested with chemotherapy drugs per ASTM D6978

Product Specifications			Physical Properties	
Gauge Thickness Measurements	MM	MIL	Before Aging	
Middle Finger	.09	3.5	Tensile Strength	42 MPa
Palm	.08	3.1	Ultimate Elongation	650%
Cuff	.07	2.8	After Aging	
Average Length	242mm (9.5")		Tensile Strength	38 MPa
	XTRA® 305mm (12.0")		Ultimate Elongation	550%

Quality Standards
Exceeds current ASTM Standards for critical defects, tensile strength, and puncture resistance.
Manufactured in accordance with Quality System ISO 9001

Reducing Space, Waste and How Much Goes to Landfill

Environmentally responsible packaging for space and waste reduction.

Environmentally Responsible Packaging		WASTE	Waste Savings Using STERLING® 9.5" Nitrile Gloves†		
9.5" Gloves	12" Gloves		Cases per Month	Saving: Lbs per Year	Saving: Tons per Year
50% More gloves per case	100% More gloves per case		25	2,200	1.0 
33% More storage space	50% More storage space		50	4,400	2.0 
39% Less waste (by weight)	59% Less waste (by weight)		100	8,600	4.0 
			150	13,000	6.0 
			200	17,400	8.0 

† Compared to KIMBERLY-CLARK® PFE 9.5" Latex Exam Gloves, waste Savings will be greater for STERLING® 12" gloves.

Best Performance

When you're in the market for a high-performance glove, turn to the proven leader. KIMBERLY-CLARK® STERLING® Nitrile Exam Gloves are a match for your demanding processes.



KIMBERLY-CLARK® STERLING® Nitrile Exam Gloves								
Description	Color	X-Small	Small	Medium	Large	X-Large	Gloves/Boxes	Total/Case
KIMBERLY-CLARK® STERLING® Nitrile Exam Gloves 9.5" Ambi	Sterling	55085	55086	55087	55088	55089	150/10 (140/10)	1500 (XL = 1400)
KIMBERLY-CLARK® STERLING® NITRILE-XTRA® Exam Gloves 12" Ambi	Sterling	53137	53138	53139	53140	53141	100/10	1000


Specific chemical resistance data can be found at: www.kc-safety.com/chemicalbarrierdata

Our Guarantee


Your total satisfaction means everything to us. If, for any reason, our products do not meet your expectations, Kimberly-Clark will reimburse you† for your initial purchase, via FREE product, for up to \$1,000. For more information on KIMBERLY-CLARK PROFESSIONAL®, visit us online at www.kcprofessional.com, contact your Kimberly-Clark Sales Representative, or call us at 1-888-346-GOKC (4652).

‡ Guarantee extended to consuming end-user accounts only.

Reduce Today, Respect Tomorrow is the KIMBERLY-CLARK PROFESSIONAL® approach to sustainability. It begins with the understanding that the way we use resources today shapes the world of tomorrow. And it has led us to focus on reducing consumption at every stage of the product lifecycle – from design and manufacture to distribution and disposal. Reduction is the key to lowering the environmental impact of our activities as well as those of customers. It has also been crucial in achieving our number one position five years in a row on the Dow Jones Sustainability Index. To learn more about Reduce Today, Respect Tomorrow and how we can reduce consumption in your business, visit www.kcpreducetoday.com



**Reduce Today
Respect Tomorrow**



Kimberly-Clark warrants that its products comply with K-C's standard specifications as of the delivery date to K-C's authorized distributors/direct purchasers. Except to the extent prohibited by applicable law: 1) **THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE;** 2) K-C is not liable for any kind of special, incidental, or consequential damages; and 3) K-C's liability for breach of contract, tort, or other cause of action shall not exceed the product purchase price. Purchasers and users are deemed to have accepted the above warranty and limitation of liability, and cannot change the terms by verbal agreement or by any writing not signed by K-C. To the extent required by applicable law, K-C does not limit its liability for death/injury resulting from K-C's negligence.