

CLEARVUE™ VLXL

COMPREHENSIVE VIDEO
LARYNGOSCOPE



CLEARVUE™ VLXL

The Infinium ClearVue™ VLXL

Video Laryngoscopes

The Infinium ClearVue™ video laryngoscope adopts state-of-the-art image technology and possesses the advantage of a 2.0 megapixel full view camera with an 8 inch high-resolution monitor. The ClearVue™ features a rechargeable li-ion battery for mobile use. The HD optics offer a unique anti-fog and high resolution image. The overall design is compact and ergonomic. A high-resolution view of the glottis will enable clinicians first-attempt success while minimizing any chance of a complication during the intubation process.

The Infinium ClearVue™
Video Laryngoscopes are:

CLEAR AND INNOVATIVE

- High Resolution 8 Inch Screen
- Mobile Stands Available
- 2 Models (Reusable Blade or Disposable Blade)
- Rigid Steel Construction of Hand Pieces
- Quick installation and Removal of Blades
- Quick Shot Camera Button for Photo and Video Recording
- 4 GB Data Storage (Approx. 80 Minutes of Stored Video)
- Use on Battery or AC Power

COMPREHENSIVE BLADE OFFERING

REUSABLE BLADE MODEL

- Macintosh Blades Size 1,2,3,4,5 Available
- Macintosh Curved “D” Blades Size 1,2,3,4,5 Available
- Miller Blades Size 0,1 available
- The ClearVue™ Comes standard with Macintosh Size 2,3, and 4 Blades



DISPOSABLE BLADE MODEL

- Macintosh Blade Size 1,2,3,4,5 Available
- Macintosh Channeled Blade Size 1,2,3,4,5 Available

DURABLE

- 2 Years Warranty
- Screen Reinforced with Anti-Shatter Protection
- Light Weight Long Lasting Rechargeable Li-ion Battery

HIGH-DEFINITION


- 8 Inch HD Display
- Anti-Fog Optics Lens
- 2.0 Mega Pixel Camera With Video

ClearVue™ Reusable Blade Model-Blades Available:



The Reusable blade model comes standard with Mac Blades sizes 2, 3, and 4

Optional Blades Available Are:



Miller O & 1 Blades



Mac 1 & Mac 5 Blades



Optional Curved Mac "D" Blades Sizes 1,2,3,4,5



ClearVue™ Disposable Blade Model-Blades Available:



The Disposable blade-model comes standard with Mac Blades sizes 1, 2, 3, 4, and 5

CLEARVUE® VLXL TECHNICAL SPECIFICATIONS:

DISPLAY

Size	8 Inch Full View
Monitor	LCD; Scratch-resistant and anti-glare screen
Resolution Ratio	1024*768
Video Refresh Rate	30FPS
Illumination	LED

LARYNGOSCOPE BLADE CAMERA INCLUDED

Resolution Ratio	2.0M Pixels
Field Angle	66°
illumination	≥800LUX

BATTERY

Type	Rechargeable Lithium Battery
Long working time	≥120minutes
Voltage	3.7V
Capacity	3200mAh
Charging time	4 hours

POWER ADAPTER

Input	100-250V,50/60Hz
-------	------------------

TRANSPORT/ STORAGE

Temperature	-10 ~+45
-------------	----------

CONDITION

Humidity	≤93
AP	500hPa to 1060hPa

WORKING CONDITION

Temperature	5 ~+40
Humidity	30 ~85
AP	700hPa to 1060hpa

ANTI-FOG FUNCTION

Auto anti-fog function upon powering on

PHOTO/VIDEO FUNCTION

One click image capture 4G internal memory card for video and images USB connection cable for uploading images and video

HANDLE

ABS material - Comfortable ergonomic short handle design

CLEARVUE ORDER INFORMATION & PART NUMBERS

000.300.0004	ClearVue VLXL, 8-in Display Laryngoscope Kit with a set of disposable blades Size 1,2,3,4,5	261.300.0013	Reusable D Blade (60) Mac 3
000.300.0009	ClearVue VLXL 8-in Display Laryngoscope Kit with a set of Reusable Blades Size 2,3,and 4	261.300.0015	Reusable D Blade (60) Mac 5
261.300.0001	Laryngoscope reusable blade No.1	306.300.0008	Laryngoscope disposable blade No.1
261.300.0002	Laryngoscope reusable blade No.2	306.300.0009	Laryngoscope disposable blade No.2
261.300.0003	Laryngoscope reusable blade No.3	306.300.0010	Laryngoscope disposable blade No.3
261.300.0004	Laryngoscope reusable blade No.4	306.300.0011	Laryngoscope disposable blade No.4
261.300.0005	Laryngoscope reusable blade No.5	306.300.0012	Laryngoscope disposable blade No.5
261.300.0008	ClearVue Reusable Blade MILLER Size 0	306.300.0013	Laryngoscope Disposable Channeled Blade No.1
261.300.0009	ClearVue Reusable Blade MILLER Size 00	306.300.0014	Laryngoscope Disposable Channeled Blade No.2
261.300.0010	ClearVue Reusable Blade MILLER Size 1	306.300.0015	Laryngoscope Disposable Channeled Blade No.3
261.300.0011	Reusable D Blade (60) Mac 1	306.300.0016	Laryngoscope Disposable Channeled Blade No.4
261.300.0012	Reusable D Blade (60) Mac 2	306.300.0017	Laryngoscope Disposable Channeled Blade No.5