



OMNI™ PRO X INFECTION CONTROL COMPLETE

M8101HV

M4L Ballistic Armored PROTECTION

**BLS/ALS Total System (TS2 Ready™)
with Optional M4L Ballistic Armored Protection**

*Refer to Specs Below

The OMNI™ PRO X is the perfect BLS/ ALS bag for the first responder who wants an organized solution for their equipment with options, and with complete infection control. Constructed internally and externally from non-absorbing wipe-clean materials the OMNI™ PRO X is specifically designed to easily clean and remove blood borne pathogens. The OMNI™ PRO X is the foundation for a system of TS2-Ready™ first responder bags that are interchangeable, enabling the first responder to attach additional, optional modules to the bag's sides depending on the team's requirements. Now available with NIJ certified Level IIIA ballistic M4L Armored Ballistic Protection, the OMNI™ PRO can withstand multiple rounds including a .44 MAG, providing protection so you can advance quicker onto an unsecure shooter scene to quickly aid and evacuate the injured.

SolidGrip, large main carry handle

Large EZ-Pull zipper pullers allow easy access

Infection control materials designed to easily clean and remove blood borne pathogens

Metal matte silver finished hardware

TS2 Ready™ hideaway zippers for optional side modules

Large front tape plate area for ID customization

Reflective paneling

Shoulder strap with metal hardware INCLUDED

Internal materials designed to easily clean and remove blood borne pathogens



Dual large external zippered pockets with mesh lined compartments

All YKK Zippers

Inner foam construction maintains shape and protects contents

Reflective logos

Laser cut MOLLE panel for modules or accessories

Front MOLLE panel for accessory customization



Bright internal materials for easy-view of supplies

Customized internal main compartment using SLIDERS™ System

Hideaway adjustable, padded, contoured yoke backpack strap harness system easily distributes weight for maximum comfort



Fits 4+ TS2 Ready™ modules internally



Adjustable sternum and waist belt, lumbar support

Zippered internal compartment with clear view pockets



Durable, lined back and bottom resists water and wear

COMFORT-GRIP carry handles on sides of bag



* M4L BALLISTIC CAPABILITIES

NIJ Standard 0108.01 Threat Protection (Level IIIA)

- 9mm Luger 124gr FMJ
- .44 MAG 240gr Lead SWC GC

Special Threat Protection

- Federal 9mm Luger 124gr Hydra-Shok JHP
- Federal .40 S&W 165gr Hydra-Shok JHP
- Federal .45 AUTO +P 200gr GDHP

Ballistic Technology Designed for MERET by Angel Armor



Product Dimensions: 15" x 22" x 9.5"
Product Weight: 11lbs 4oz standard,
M4L Lifetime Warranty

Available TS2 Ready™ Modules

M8101A-F

AIRWAY PRO X
Intubation
Tri-Fold Module



M8101C-F

XTRA FILL PRO X
Multi-Use Module



M8101J-F

MEDKIT PRO X
Medications Module



OMNI PRO X also available in:

Blue
Infection Control
COMPLETE



M8101NB

Red
Infection Control
COMPLETE



M8101F

Tactical Black
Infection Control
COMPLETE



M8101TB

877.222.0200 | Ph: 714.259.4700 | Fx: 714.259.4749 | www.MeretUSA.com | info@meretusa.com

BUILT FOR THE PROFESSIONAL

MEDICAL * EMERGENCY * RESPONSE * EQUIPMENT & TECHNOLOGY