



FEATURES

Accommodate a large variety of sharps sizes.

Temporary and final closure.

Built-in handles for easy transport (excluding phlebotomy).

Translucent containers and lids help identify fill level.

Nestable to save space and reduce shipping costs.

Horizontal drop on most units maximizes container volume.

Can be autoclaved or incinerated for final disposal.

Multilingual labeling ensures user safety.

PHLEBOTOMY & MULTI-USE SHARPS CONTAINERS

Designed for use in restricted areas such as nurses' stations, med carts, laboratories, ORs and pharmacies. Can be wall-mounted or used freestanding on tabletop or cart.





100 030
100 040

1-QUART PHLEBOTOMY CONTAINER

100/CASE

Translucent Red

Translucent Yellow

- Fits phlebotomy trays
- Dual openings; vertical drop
- Dimensions: 6 ⁷/₁₆" H x 4 ¹/₄" L x 4 ¹/₄" W
- *Uses holder #400 to secure container on tray or tabletop*



400 000

PHLEBOTOMY HOLDER

20/CASE

- Translucent White allows visibility of container label
- Adhesive on bottom secures container to tabletop or tray
- Dimensions: 6 ⁵/₁₆" H x 4 ¹/₄" L x 4 ¹/₄" W
- *For use with phlebotomy container #100*



102 020
102 030
102 040

2-GALLON SHARPS CONTAINER

30/CASE

Translucent Beige

Translucent Red

Translucent Yellow

- Accommodates syringes up to 60cc
- Built-in handles for easy transport
- Dimensions: 8 ⁵/₈" H x 11 ⁵/₈" L x 7 ³/₄" W
- *Uses bracket/key #435 and 440 for wall mounting*



303 020
303 030
303 040

3-GALLON SHARPS CONTAINER

12/CASE

Translucent Beige

Translucent Red

Translucent Yellow

- Large capacity accommodates syringes up to 60cc
- Built-in handles for easy transport
- Dimensions: 13 ¹/₂" H x 13 ⁷/₈" L x 6 ⁷/₈" W
- *Uses bracket #423 or 425 for wall mounting*



FEATURES

Rotating lid design prevents access to container contents.

Horizontal drop allows maximum capacity of sharps.

Translucent containers and lids help identify fill level.

Containers lock for final disposal.

Accommodate syringes up to 60cc.

Nestable to save space and reduce shipping costs.

Can be autoclaved or incinerated for final disposal.

Multilingual labeling ensures user safety.

PATIENT ROOM & EXAM ROOM SHARPS CONTAINERS

Designed for use in patient rooms, physicians' offices, clinics and other high-traffic areas.





-  125 020
-  125 030
-  125 040

3-QUART WALLSAFE® SHARPS CONTAINER 24/CASE

- Translucent Beige**
- Translucent Red**
- Translucent Yellow**

- Rotating cylinder top allows one-handed disposal
- Ideal for low-volume areas
- Dimensions: 8" H x 11" L x 4 1/4" W
- *Requires use of bracket/key #410 and 415 for wall mounting*



-  150 020
-  150 030
-  150 040

5-QUART WALLSAFE® SHARPS CONTAINER 24/CASE

- Translucent Beige**
- Translucent Red**
- Translucent Yellow**

- Rotating cylinder top allows one-handed disposal
- Dimensions: 12" H x 11" L x 4 1/4" W
- *Requires use of bracket/key #410 and 415 for wall mounting*



-  175 020
-  175 030
-  175 040

5-QUART SHARPS CONTAINER 32/CASE

- Translucent Beige**
- Translucent Red**
- Translucent Yellow**

- Rotating cylinder top allows hands-free disposal
- Built-in handles for easy transport
- Dimensions: 10" H x 11" L x 5 1/4" W
- *Uses bracket/key #435 and 440 or cabinet 1051 for wall mounting*



-  333 020
-  333 030
-  333 040

3-GALLON SHARPS CONTAINER 12/CASE

- Translucent Beige**
- Translucent Red**
- Translucent Yellow**

- Rotating cylinder top allows hands-free disposal
- Built-in handles for easy transport
- Dimensions: 15" H x 13 7/8" L x 6 7/8" W
- *Uses bracket #423 or 425, or cabinet 1031 for wall mounting*



SHARPS WALL CABINET SYSTEM

Secures sharps containers in a sturdy, locked enclosure. Made of tough ABS plastic for durability and ease of cleaning. Neutral color blends with any décor. Tinted window allows for easy viewing of fill level. Available in two sizes: 3-Gallon and 5-Quart. Optional glove box holder keeps gloves readily available.





1051

5-QUART WALL CABINET

1/CASE

- Two universal keys provided
- Cabinet back accommodates most competitive drilling/wall-mounting patterns
- Dimensions: 11 1/2" H x 13" L x 6" W
- *For use with sharps container #175*



1031

3-GALLON WALL CABINET

1/CASE

- Two universal keys provided
- Cabinet back accommodates most competitive drilling/wall-mounting patterns
- Dimensions: 16 3/4" H x 15" L x 7" W
- *For use with sharps container #333*



1041

GLOVE BOX HOLDER

1/CASE

- Snaps easily onto bottom of cabinets
- Dimensions: 5 1/2" H x 11 5/16" L x 4 1/8" W
- *For use with wall cabinets #1051 and 1031*



1060

SHELF FOR 3-GALLON CABINET

12/CASE

- Dimensions: 13 7/8" L x 4 1/2" W
- *Converts cabinet #1031 to accommodate 5-quart sharps container #175*





FEATURES

Compact design uses minimum floor space.

Large opening lids accommodate a variety of sharps.

Lids provide temporary and final closure.

Nestable to save space and reduce shipping costs.

Molded-in handles for easy lifting and transport.

Multilingual labeling ensures user safety.

LARGE VOLUME SHARPS CONTAINERS

Designed for high-volume use in operating rooms, pharmacies, laboratories, emergency rooms, dialysis rooms, and other areas. A rolling cart is available for convenient mobility and transport.





■ 108 030
■ 108 040

8-GALLON SHARPS CONTAINER

10/CASE

Red

Yellow

- Separate round opening accepts small items
- Molded-in handles for easy lifting and transport
- Hinged lid accommodates large items
- Dimensions: 15 7/8" H x 16 1/2" L x 11 13/16" W
- *Uses cart #430 for mobility*



■ 111 030
■ 111 040

11-GALLON SHARPS CONTAINER

6/CASE

Red

Yellow

- Separate round opening accepts small items
- Molded-in handles for easy lifting and transport
- Hinged lid accommodates large items
- Dimensions: 22 1/2" H x 16 1/2" L x 11 13/16" W
- *Uses cart #430 for mobility*



FEATURES

Gasketed lids and absorbent spill pads provide leak resistance.

Available in three sizes.

Clearly labeled with chemotherapy warnings.

Nestable to save space and reduce shipping costs.

Multilingual labeling ensures user safety.

Note: Chemotherapy containers must be incinerated.

CHEMOTHERAPY SHARPS CONTAINERS

Designed for safe collection and containment of trace chemotherapy waste including needles, syringes, IVs, empty or partially empty vials, as well as gloves and gowns. The EPA defines "trace" amounts of waste as containing no more than 3% by weight of the original contents.





202 004
 202 040

2-GALLON CHEMOTHERAPY CONTAINER 30/CASE

White w/gasket lid/absorbent pad
Yellow w/gasket lid/absorbent pad

- Screw-on cap; Posts on lid for holding cap in place when in use
- For use on tabletop or wall-mounted
- Dimensions: 9" H x 11 5/8" L x 7 3/4" W
- *Uses bracket/key #435 and 440 for wall mounting*



208 004
 208 040

8-GALLON CHEMOTHERAPY CONTAINER 10/CASE

White w/gasket lid/absorbent pads
Yellow w/gasket lid/absorbent pads

- Separate round opening accepts small items
- Hinged lid accommodates large items
- Molded-in handles for easy lifting and transport
- Dimensions: 15 7/8" H x 16 1/2" L x 11 13/16" W
- *Uses cart #430 for mobility*



211 004
 211 040

11-GALLON CHEMOTHERAPY CONTAINER 6/CASE

White w/gasket lid/absorbent pads
Yellow w/gasket lid/absorbent pads

- Separate round opening accepts small items
- Hinged lid accommodates large items
- Molded-in handles for easy lifting and transport
- Dimensions: 22 1/2" H x 16 1/2" L x 11 13/16" W
- *Uses cart #430 for mobility*



FEATURES

Gasketed lids and absorbent pads provide leak resistance.

Available in three sizes.

Nestable to save space and reduce shipping costs.

NOT for sharps.

Multilingual labeling ensures user safety.

Note: Pharmaceutical waste must be incinerated.

PHARMACEUTICAL WASTE CONTAINERS

Designed for the disposal of expired or unused non-hazardous pharmaceuticals such as antibiotics, IVs, ibuprofen, etc. Ideal for use in pharmacies, nursing stations and med carts. Blue color and special labeling easily distinguishes from other disposal containers for sharps, chemotherapy and RCRA waste.





4002 050

2-GALLON PHARMACEUTICAL CONTAINER

30/CASE

Blue w/gasket lid/absorbent pad

- Screw-on cap; Posts on lid for holding cap in place when in use
- For use on tabletop or wall-mounted
- Dimensions: 9" H x 11 ⁵/₈" L x 7 ³/₄" W
- *Uses bracket/key #435 and 440 for wall mounting*



4008 050

8-GALLON PHARMACEUTICAL CONTAINER

10/CASE

Blue w/gasket lid/absorbent pads

- Separate round opening accepts small items
- Hinged lid accommodates large items
- Molded-in handles for easy lifting and transport
- Dimensions: 15 ⁷/₈" H x 16 ¹/₂" L x 11 ¹³/₁₆" W
- *Uses cart #430 for mobility*



4011 050

11-GALLON PHARMACEUTICAL CONTAINER

6/CASE

Blue w/gasket lid/absorbent pads

- Separate round opening accepts small items
- Hinged lid accommodates large items
- Molded-in handles for easy lifting and transport
- Dimensions: 22 ¹/₂" H x 16 ¹/₂" L x 11 ¹³/₁₆" W
- *Uses cart #430 for mobility*



FEATURES

Gasketed lids and absorbent pads provide leak resistance.

Available in three sizes.

Nestable to save space and reduce shipping costs.

NOT for sharps.

Multilingual labeling ensures user safety.

Note: RCRA waste must be incinerated by a federally permitted hazardous waste incinerator.

RCRA HAZARDOUS WASTE CONTAINERS

The Resource Conservation and Recovery Act (RCRA) was enacted by the EPA in 1976 to regulate the disposal of solid hazardous chemical wastes. Bemis RCRA Hazardous Waste Containers were designed to meet the EPA requirements for disposal of P-, U- and D-listed chemicals. Designed for use in pharmacies and nursing stations, these containers are black and contain special RCRA labeling to distinguish from other waste containers.





5002 070

2-GALLON RCRA CONTAINER

30/CASE

Black w/gasket lid/absorbent pad

- Screw-on cap; Posts on lid for holding cap in place when in use
- For use on tabletop or wall-mounted
- Dimensions: 9" H x 11 ⁵/₈" L x 7 ³/₄" W
- *Uses bracket/key #435 and 440 for wall mounting*



5008 070

8-GALLON RCRA CONTAINER

10/CASE

Black w/gasket lid/absorbent pads

- Separate round opening accepts small items
- Hinged lid accommodates large items
- Molded-in handles for easy lifting and transport
- Dimensions: 15 ⁷/₈" H x 16 ¹/₂" L x 11 ¹³/₁₆" W
- *Uses cart #430 for mobility*



5011 070

11-GALLON RCRA CONTAINER

6/CASE

Black w/gasket lid/absorbent pads

- Separate round opening accepts small items
- Hinged lid accommodates large items
- Molded-in handles for easy lifting and transport
- Dimensions: 22 ¹/₂" H x 16 ¹/₂" L x 11 ¹³/₁₆" W
- *Uses cart #430 for mobility*



BRACKETS & ACCESSORIES




400 000
PHLEBOTOMY HOLDER
20/CASE


- Translucent White allows visibility of container label
- Adhesive on bottom secures container to tabletop or tray
- Dimensions: 6 5/16" H x 4 1/4" L x 4 1/4" W
- *For use with phlebotomy container #100*


410 020
MOUNTING BRACKET FOR WALLSAFE® SHARPS CONTAINERS
10/CASE


- Holds WallSafes #125 and 150 securely to wall
- Concealed by container when in use
- Mounting hardware included
- Dimensions: 6" H x 9" L x 1/2" W
- *Requires key #415 to remove container from bracket*


415 020
KEY FOR #410 MOUNTING BRACKET
10/CASE

- Unique plastic key removes WallSafes from bracket
- Dimensions: 3 1/4" H x 3/4" L x 3/4" W
- *For use with mounting bracket #410*


423 000
NON-LOCKING WIRE BRACKET FOR 3-GALLON
5/CASE

- Holds sharps containers #303 and 333 securely to wall
- Mounting hardware included
- Dimensions: 7 1/2" H x 12" L x 6" W
- *Can be mounted separately or slides onto mounting bracket #410*


425 000
LOCKING WIRE BRACKET FOR 3-GALLON
5/CASE

- Holds sharps containers #303 and 333 securely to wall
- Same as bracket #423 but with swing arm that fits over container and secures with padlock
- Mounting hardware included
- Dimensions: 9 1/2" H x 12" L x 6" W
- *Can be mounted separately or slides onto mounting bracket #410*
- *Requires padlock/key #426*


426 000
PADLOCK & KEY
1/CASE

- Dimensions: 1 3/4" H x 1 1/8" L x 1/2" W
- *For use with locking wire bracket #425*


430 030
ROLLING CART FOR 8- & 11-GALLON CONTAINERS
1/CASE

- Durable casters
- Sturdy construction
- Dimensions: 16" H x 17" L x 13" W
- *Accommodates all 8- and 11-gallon containers*


435 020
MOUNTING BRACKET FOR 5-QUART & 2-GALLON
10/CASE

- Dimensions: 6" H x 10 1/4" L x 2 1/4" W
- Mounting hardware included
- *Holds sharps container #175 5-quart and all 2-gallon containers securely to wall*
- *Requires key #440 to remove container from bracket*


440 020
KEY FOR #435 MOUNTING BRACKET
10/CASE

- Unique plastic key removes #175 5-quart and all 2-gallon containers from bracket
- Dimensions: 3 1/8" H x 1" L x 5/8" W
- *For use with mounting bracket #435*