

Med System III

Med System III set Half Set 1 SmartSite needle-free valve(s) 9" (below pumping segment) from 2-piece Male Luer Lock. L: 85 in L: 216 cm PV: 6 mL Fluid Path Sterile

Summary

SKU	28117E
Product Name	Med System III
Product Type	



Weight (lbs)	6.6000	Country of Manufacture	Mexico	Product Sub Type	No
Filtration Surface Area	No	Anti-Siphon Valve	No	Approximate Length (cm)	216
Approximate Length (in)	85	Approximate Priming Volume (mL)	6	Bag Access	No
Ball Valve in Drip Chamber	No	Blood Set	No	Not made with BPA/Polycarbonate	No
Burette	No	Case Dimensions	23.375in x 14.125in x 6.125in 6.6lb	Clamp	No
DEHP is not part of the material formulation	No	Drop	N/A Half Set	Flow Regulator	No
Hanger	No	Inner Diameter (I.D.)	No	Natural Rubber Latex is not part of the material formulation	Yes
Low Sorbing	No	Luer Lock	2-piece Male Luer Lock	Manifold	None
Number of Back Check Valves	0	Number of Connectors	0	Number of Fuses	1
Pressure Rated	No	Pressure Sensing Disc	No	Sterility	Fluid Path Sterile
Stopcock	None	Texium	No	Color of Tubing	None
T-Connector Type	None	Number of y-sites	1	Y-Site Port Type	SmartSite
Type of Connector	None	Filter	None	Tubing Size	Standardbore
CE Marking	No	PVC Free	No	Quantity	50
Unit Measure	Per case				