

Microflex®
SAFEGRIP®
 Powder-Free Latex
 Extended Cuff Exam Gloves

Item #

Small SG-375-S
 Medium..... SG-375-M
 Large..... SG-375-L
 Extra Large..... SG-375-XL



Recommended for use in the following industries:

Automotive / EMS / Healthcare / Industrial Safety / Laboratory

For the highest level of protection and durability in latex exam gloves for heavy-duty or high-risk environments, rely on **SafeGrip®** gloves.

SafeGrip® gloves are up to 4 times thicker than standard latex exam gloves for defense against many hazards, and is ideal protection for your toughest tasks.

SafeGrip® gloves are the best choice in a latex for heavy-duty/high-risk applications.

Glove meets or exceeds the emergency medical glove requirements of NFPA 1999, Standard on Protective Clothing for Emergency Medical Operations, 2008 edition.*

Recommended for use when:

- Added protection is required for extreme environments
- Extra coverage over the wrist and forearm is needed for enhanced safety and security
- Working in situations requiring intense pressure and strain on the hands

*Certification on file.

Microflex® SafeGrip® Profile			For current specification information, please visit www.microflex.com .	
	ASTM D 3578	SafeGrip®	Material: Latex	Type: Non-sterile
Length (mm/in):	230 / 9.1	300 / 11.8	Size: S - XL	Color: Blue
Cuff Thickness (mils):	n/a	7.9	Shape: Ambidextrous	Usage: For single use only
Palm Thickness (mils):	min. 2.0	11.0	Cuff: Extended Cuff length (12"), Beaded	
Finger Thickness (mils):	min. 2.0	14.2	Glove Exterior: Textured Finger Tips	
Tensile Strength:			Glove Interior: Powder-Free, Chlorinated	
Before Aging (MPa):	min. 18	34	Packaging: 500 gloves/case, by weight; 10 boxes/case	
After Aging (MPa):	min. 14	28		
Elasticity/Elongation:			50 gloves per dispenser box, by weight	
Before Aging (%):	min. 650	750		
After Aging (%):	min. 500	600		
			Caution: This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.	

The reason is the quality... the difference is Microflex®.

See for yourself why Microflex® gloves are the first choice of exam glove users nationwide. For **FREE SAMPLES**, please visit www.microflex.com. For the name of a distributor in your area, call: **(800) 876-6866**.

