

# ONLY ONE

## FL<sub>OW</sub>SAFE II EZ® DISPOSABLE CPAP SYSTEM

### CPAP System with Integrated Nebulizer

A major leap in product innovation, it takes Emergency Care to a whole new level.

- CPAP & Nebulizer System in One
- **ONLY ONE** O<sub>2</sub> Source
- Less Parts – Easy Set-up



On/Off nebulizer switch incorporated in new integrated design



[www.mercurymed.com](http://www.mercurymed.com)

# ONLY ONE

## FLOWSAFE II EZ<sup>®</sup> DISPOSABLE CPAP SYSTEM

When lives hang in the balance, you need a disposable CPAP system that lives up to the challenge.

The Flow-Safe II EZ<sup>®</sup> Disposable CPAP System from Mercury Medical<sup>®</sup> is the only system that has so many advantages, it gives you an advantage in life-saving emergencies.

In short, the Flow-Safe II EZ represents nothing short of a revolution in emergency care.



On/Off nebulizer switch incorporated in new integrated design



# SAFE. SURE. SUPERIOR. CPAP.





Clinicians can *now* deliver aerosol & CPAP therapy with just **ONE** integrated oxygen source.

### ADVANTAGE 1:

#### INTEGRATED NEBULIZER

The Flow-Safe II EZ offers the capability of an in-line nebulizer.

- Uses only one oxygen supply source
- Easy set-up, less parts

### ADVANTAGE 2:

CPAP SYSTEM CONSUMES 50% LESS OXYGEN\* THAN ORIGINAL FLOW-SAFE®

The Flow-Safe II EZ conserves oxygen while maintaining high FiO<sub>2</sub> delivery. Saves resources and represents a major advantage on long transports. Increasing flow may be necessary when activating the nebulizer.



### ADVANTAGE 3:

#### BUILT-IN MANOMETER & PRESSURE RELIEF VALVE

Only the Flow-Safe II EZ system features a built-in manometer for verified pressure readings. No assembly of separate apparatus. And the pressure relief valve automatically adjusts to avoid excess pressure.

### ADVANTAGE 4:

#### ADVANCED MASK DESIGN

The Flow-Safe II EZ deluxe face mask is lighter, easier to handle and is designed to form a better anatomical seal. The elastic head harness is easy-to-place with velcro straps that easily adjust for patient comfort.



### ADVANTAGE 5:

ACCEPTS STANDARD CO<sub>2</sub> SAMPLING PORT  
CO<sub>2</sub> Sampling Port Fits In-Line (supplied separately)

PRSRT STD  
U.S. POSTAGE  
PAID  
Permit #579  
St. Petersburg, FL

Clinicians can *now* deliver aerosol & CPAP therapy with just **ONE** integrated oxygen source.



**FLOWSAFE II EZ**  
DISPOSABLE CPAP SYSTEM

### ORDERING INFORMATION

Part #	Description	Packaged
#10-57318	Flow-Safe II EZ CPAP System, large adult deluxe full face mask with ports and EZflow Max Nebulizer	5/Box
#10-57319	Flow-Safe II EZ CPAP System, medium/small adult deluxe full face mask with ports and EZflow Max Nebulizer	5/Box
#10-57320	Flow-Safe II EZ CPAP System, small/child deluxe full face mask with ports and EZflow Max Nebulizer	5/Box
#10-57321	Flow-Safe II EZ CPAP System, without mask with EZflow Max Nebulizer	5/Box

***A Revolution in Emergency Care***

**MercuryMedical**<sup>®</sup>

11300 - 49th Street North • Clearwater, Florida 33762-4807 USA

**ORDERING INFORMATION: 800.835.6633 or 727.573.0088**

**Fax: 727-573-6040**

**[www.mercurymed.com](http://www.mercurymed.com)**

Patents: US 8,522,618 B1 / US 6,338,443 B1 /  
US 9,370,635 B2 / US 9,511,202 / Other patents pending