

FilterLine® Selection Guide

Sampling solutions for Microstream® enabled monitors

NON-INTUBATED PATIENTS

Long term use

Smart CapnoLine H Plus



Key Applications
Medical surgical floor,
Critical Care, Sleep Lab

Oral-Nasal

Smart CapnoLine H



Key Applications
Medical surgical floor,
Critical Care

Oral-Nasal

CapnoLine H



Key Applications
When nasal sampling is
preferred, Critical Care,
Sleep Lab

Nasal

Short term use

Smart CapnoLine® Plus



Key Applications
Procedural Sedation, Lower
GI procedures, OR, MAC,
EMS, ED, Rapid Response
Team

Oral-Nasal

Smart CapnoLine



Key Applications
Procedural Sedation, EMS,
ED, Rapid Response Team

Oral-Nasal

Smart CapnoLine Guardian™



Key Applications
EGD, ERCP,
Bronchoscopy,
Transesophageal
Echocardiogram

Oral-Nasal

O₂/CO₂ Nasal FilterLine



Key Applications
When nasal sampling is
preferred

Nasal

Nasal/NIV Line™



Key Applications
EMS, ED

Nasal

INTUBATED PATIENTS

Long term use

FilterLine H Set




Key Applications
Critical Care, Humidified
environments

VitaLine™ H Set



Key Applications
Critical Care, High
ambient humidity

FilterLine Set for Infant/Neonates



Key Applications
Critical Care, Humidified
environments

VitaLine H Set for Infant/Neonates



Key Applications
Critical care, Humidity-
controlled incubators

Short term use

FilterLine® Set



Key Applications
OR, EMS, ED, Rapid
Response Team, Transport



Oridion®

www.oridion.com

Intubated patient sampling lines



FilterLine® Set

| | |
|---------------------------|---------|
| Adult/Pediatric | XS04620 |
| Adult/Pediatric 100 count | 010579 |
| Adult/Pediatric Long | 007768 |





FilterLine® H Set

| | |
|---------------------------|---------|
| Adult/Pediatric | XS04624 |
| Adult/Pediatric 100 count | 010580 |
| Adult/Pediatric Long | 007737 |
| Infant/Neonates | 006324 |
| Infant/Neonates Long | 007738 |



VitaLine™ H Set

| | |
|-----------------|--------|
| Adult/Pediatric | 010787 |
| Infant/Neonates | 010807 |



Product features & configurations

| Short term duration | Long term duration | Adult size | Pediatric size | Infant/Neonate size | 25 unit packaging | 100 unit bulk packaging | Long length (15 ft.) |
|---------------------|--------------------|------------|----------------|---------------------|-------------------|-------------------------|----------------------|
| duration | | size | | | packaging | | |
| • | | • | • | | • | | |
| • | | • | • | | | • | |
| • | | • | • | | • | | • |
| | • | | | | | | |
| | • | • | • | | • | | |
| | • | • | • | | | • | |
| | • | • | • | • | • | | • |
| | • | • | • | | • | | • |
| | • | | | | | | |
| | • | • | • | | • | | |
| | • | | | | | | |

All FilterLine® breath sampling products are for single patient use only.

All FilterLine breath sampling products are latex free.

Packaged in 25 units unless otherwise specified.

Standard FilterLine Length: **200 cm, 6.5 ft.** Long FilterLine Length: **400 cm, 13 ft.**

Adult/Pediatric products are intended for use with ETT tube size > 4.5 mm.
The adapter dead space is < 6.6 cc and the weight is 7.0 g.

Infant/neonate products are intended for use with ETT tube size ≤ 4.5 mm.
The adapter dead space is < 0.5 cc and the weight is 3.8 g.

FilterLine breath sampling products, with the exception of infant/neonate intubated patient sampling lines (Infant/Neonate FilterLine H part nos 006324 / 007738 and VitaLine H set part no. 010807), are compatible in the MRI environment.

Specifications are subject to change without notice.

Microstream, FilterLine, CapnoLine, and Capnostream are registered trademarks and NIV Line, Smart CapnoLine Guardian, and VitaLine are trademarks of Oridion Medical 1987 Ltd. Oridion products are covered by one or more of the following US patents: 6,422,240; 5,857,461; 6,437,316; 6,926,005; 5,657,750; 7,383,839 and their foreign equivalents. Additional patent applications pending.



Oridion Capnography Inc.
160 Gould Street
Suite 205
Needham MA 02494
USA
Tel: (888) ORIDION
Fax: (781) 453-2722
www.oridion.com

Oridion Medical 1987 Ltd.
7 Hamarpe Street
97774
Jerusalem
Israel
Tel: +972 2 589 9158
Fax: +972 2 586 6680